

Edge Mode+ ENHANCED



Now with Automatic Adaptation

The industry's first-of-its-kind intelligent, on-demand noise management and comfort feature will now continue to scan, adapt, and optimize as the wearer changes environments or as the acoustic environment changes — ensuring the most appropriate settings are always applied in a dynamic world and eliminating the need for multiple programs.

Automatic adaptations are made in any Edge Mode+ setting - Best sound, Enhance speech, or Reduce noise.

GENESIS^{AI} New and previously fit 2.4 GHz devices with a firmware update.





Figure 1: Edge Mode+ control options in the My Starkey mobile app

Edge Mode+ has been shown to improve speech recognition and decrease listening effort.¹



Where to find it

Genesis AI		
My Starkey mobile app	User Control	Remote Control

Patient Feedback



Edge Mode+ is impressive. It worked so good. It was almost unsettling at first because I felt like I wasn't hearing well, and then what I realized is I wasn't hearing the noise well, but I was understanding voices perfectly.





Edge Mode+ Set-up & Use Instructions

Edge Mode+ is available in any wireless Genesis AI product. It is defaulted as an option on the homescreen of the My Starkey app, but may also be configured in the hearing aid User Controls and/or the Remote Control 2.0 Favorite Button. See below for details on engaging Edge Mode+.

Hearing aid user control

In Pro Fit select **User Controls** on the Flyout Menu .


Drag  control from the left and drop to either **Short Press** or **Double Tap** gesture on the right.



Remote control using the Favorite Button

In Pro Fit select **Accessories** on the Flyout Menu .

Use the dropdown menu below the Remote image to set the Favorite button configuration.



My Starkey mobile app

1. Open the My Starkey mobile application.
2. Tap the **Edge Mode+** icon on the home screen.
3. Best Sound is automatically applied; for further personalization based on the patient's listening intent: speech or reduced noise, the patient can select either Enhanced Speech or Reduce Noise.
4. Tap **Save**.
5. To exit Edge Mode+, select **Programs** and select the desired program.

