

**QUICKTIP**

# Hear Clear Receiver Wax Guards

## Product Availability: RIC and Custom Devices

The Hear Clear earwax protection system uses disposable wax guards. The innovative wax guards prevent earwax accumulation in the hearing device receiver. Hear Clear application sticks are yellow.

Routine maintenance or replacement of Hear Clear Receiver Wax Guards helps to ensure optimal performance of hearing devices.


**NOTE:** *Hear Clear Receiver Wax Guards are not compatible with other manufacturers' devices.*

## Insertion and Removal


- 1 Insert the empty end of the yellow application stick into the used wax guard in the hearing device.

**NOTE:** *Push the yellow stick directly into the used Hear Clear. There is no need to twist the yellow stick.*

- 2 Pull outward on the stick to remove the used wax guard.
- 3 Use the opposite end of the yellow stick to firmly insert the clean wax guard into the hearing device.
- 4 Pull outward to remove the yellow stick and discard.



RIC Devices



Custom Devices

## Use

- Hear Clear Receiver Wax Guards are available to order in packs of 8 and 48.
- Packs of 8 will ship with all RIC products, custom products, RIC receivers, and RIC Absolute Power (AP) or Custom Cased receivers.
- Orders include: Blue packs of 8 yellow wax guards and instructions on how to insert and replace them.