### T1130-01-03

# Evaluating the Versatility of a New High Molecular Weight HPC in Dual Active Diabetic Controlled Release Dosages

### **CONTACT INFORMATION: Teslin.botoy**@ashland.com

# PURPOSE


Dual-active controlled-release formulations for Type Il diabetes provide several benefits, including reduced dosing frequency, improved glycemic control, and improved patient compliance. Common hydrophilic matrix-forming polymers such as Hydroxypropyl Methylcellulose (HPMC) usually comprise 20-30% of the weight of a tablet to achieve an acceptable release profile and require intra-granulating binders to prevent the burst effect when dealing with highly water-soluble drugs such as Metformin. This can result in large, hard-to-swallow tablets with complex processing.

## **OBJECTIVES**

This study evaluates the process versatility and the tablet size reduction capabilities of the new high molecular weight hydroxypropyl cellulose, klucel xtend<sup>™</sup> to match the release profiles of both, a bilayer dosage and a metformin HCl core tablet coated with an immediate-release drug substance for dual diabetic drug delivery.


### **METHODS**

The polymer gel strength of the new hydroxypropyl cellulose (HPC), klucel xtend<sup>™</sup> and hydroxypropyl methylcellulose (HPMC) K100M was determined using a Discovery Hybrid Rheometer (DHR-3) with a Peltier heating stage. Dynamic mechanical analysis was used to test the thermoplasticity of both polymers.





A Instruments iscoverv Hybrid Rheomet Peltier heating stage 25 °C, linear rate of 10  $\mu$ m/s for 500 s



Probe goes down (negative direction) and measures

resistance force The higher the resistance, the stiffer is the gel laver

strain = distance / original height stress = (resistance) force / area stress

 $\frac{1}{1}$  = compression modulus strain

Metformin HCl was wet granulated with HPC at a w/w ratio of 85:15. dried in an oven at 65°C overnight, then milled with a 45mesh Quadro<sup>®</sup> Comil<sup>®</sup> screen. Granules were then mixed with extra-granular materials and compressed as either a monolayer core or a bi-layer tablet with the drug dapagliflozin in the IR layer. The monolayer cores were later coated with a barrier layer, followed by a drug layer of empagliflozin using a Freund Vector LDCS<sup>®</sup> Tablet Coater. Both tablets were then cosmetically filmcoated.

Magnes

Aquarius™ Total tabl

All tablets, including the bi-layer Xigduo<sup>®</sup> XR and core/coated dual drug Synjardy<sup>®</sup> XR were characterized for dimensions, hardness, friability, and dissolution release profiles. Dissolution of each formulation was conducted in 900 mL of pH 6.8 0.05M phosphate buffer maintained at 37°C with USP Apparatus 1 (baskets) at 100 RPM. Samples were taken at intervals of up to 24 hours; quantitation was by HPLC with UV detection for all three drugs. ered trademark of Astra Zeneca US, Synjardy<sup>®</sup> XR is a registered trademark of Boehringer Ingelheim and Glucophage<sup>®</sup> is a registered idemark of Merck Santé S.A.S., an associate of Merck kGaA of Darmstadt, Germany, Licensed to Bristol-Myer Squibb Company

Teslin Botoy, Brian Huebner, Quyen Schwing, Kapish Karan, and Thomas Dürig

Ashland Specialty Ingredients, Wilmington, DE 19808

и.					
gredient	Tablet Formulation		Ingradiant	<b>Tablet Formulation</b>	
	W/W (%)	Wt (mg)	Ingredient	W/W (%)	Wt (mg)
Immediate release layer			Controlled release layer		
oagliflozin	10	10	Metformin	84	1040
el EXF hpc	3	3	Klucel xtend hpc	15	186
asdone XL-10 spovidone	10	10	Magnesium stearate	1.0	12
с С РН 102	76	76	Total	100	1238
sium stearate	1	1	Separating layer		
Total	100	100	HPMC E5	2	25
Controlled release layer			Immediate release drug coating		
etformin	84	1029	HPMC E5	1.14	14
el xtend hpc	15	184	empagliflozin	0.86	1
sium stearate	1	12	Cosmetic Co		
Total	100	1226			
PVA Film Coating (mg)	3	40	Aquarius™ Preferred HSC Film Coating (mg)	3	37
et Weight (mg) 1365		Total tablet Weight (mg) 1325			
ended for 10 minutes in a turbula blender ble			batch of metformin was blended with Klucel™ xtend		
-					
of metformin was blended with Klucel™ xtend		Metformin blend wet granulated at 750 rpm 1500 rpm chopper			
-					
nin blend wet granulated at 750 rpm 1500 rpm chopper			each granulation was dried at 65°C milled		
granulation was d milled	anulation was dried at 65°C milled comp		oression CR blend was blended with magnesium stearate for 2 minutes		
, both IR and CR blends were d with magnesium stearate and		CR compressed into tablets			
blended for 2 mir	nutes			+	
compressed in bilayer tablets, ed using Aquarius™PVA and tested for dissolution			3 layers of coating, first HPMC E5, then HPMC E5 and empagliflozin and finally coated using Aquarius prefer HSC and tested for dissolution		

integrity.

tablet volume reduction by -20.43%

