



AUDIBEL

Product Portfolio Guide

Summer 2026



AUDIBEL

Tech that's yesterday's thinking **or future driven?**

From enhancements to our in-demand Aris AI to a new essential line featuring two new BTEs, Audibel is setting the speed of innovation with technology that lets you fit AI-powered performance on every patient—all using the same fitting software and app.

ASK MORE from your technology partner

ENHANCED

ARIS AI

Tech Tier/Channels/Bands

24 / 20 / 16

AURACAST

CROS COMPATIBLE



mRIC R



RIC RT



RIC 312



CIC



ITC R



ITE R



- Waterproof*
- Rechargeable
- LED Status Indicator Light
- Wireless
- Telecoil

*Up to 1m.

Available Colors

- Standard Colors**
- White
 - Tech Black
 - Graphite Gray
 - Silver
 - Chestnut
 - Caramel
 - Beige

- Custom Colors**
- Faceplate Colors**
- Pink
 - Light Brown
 - Medium Brown
 - Chestnut
 - Dark Brown
 - Black
- Shell Colors**
- Clear
 - Pink
 - Light Brown
 - Blue/Red
 - Black

SnapFit 2.0 Stock and Custom Cased Receivers



Comfort Bud Open
5 mm, 6 mm, 7 mm,
9 mm, 11 mm



Comfort Bud Occluded
4 mm, 5 mm, 6 mm,
7 mm, 9 mm, 11 mm



Power Dome
8-10 mm, 10-12 mm,
12-14 mm



Custom Cased Receiver
Custom shell with
embedded receiver



**Stock Silicone
Helix RIC Earmold**






**Stock Hollow
Silicone RIC Earmold**



**Stock Silicone
RIC Earmold**

Fitting Ranges



Visit AudibelPro.com/AriseAI for full technical data sheets.



Engineered with our industry-leading AI, our G3 Gen AI Neuro Processor, and the world's first DNN-powered directionality system, DNN 360, Aris AI is Audibel's most advanced hearing technology yet. It's designed to come closer to restoring the normal sense of hearing than any hearing technology before it.

For more Aris AI details,
visit AudibelPro.com/AriseAI



ARIS AI

Tech Tiers | Channels | Bands

24

20

16

Edge Mode+				
Automatic & User Control Personalization App Control		●	○	○
Speech Optimization				
Consonant Brightness Voice AI*		●	○	○
Sound Manager				
DNN 360	DNN Directionality**	●	○	○
	DNN Spatial Awareness**	●	○	○
	Quiet	Maximum expansion	Moderate expansion	Mild expansion
	Machine Noise Reduction	Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction
	Speech in Noise Reduction	Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction
	Transient Noise Reduction	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 6 dB of reduction
	Wind Reduction	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction
Additional Key Features				
Gen AI Smart Assistant & TeleHear AI		●	○	○
Translate		●	○	○
ALL NEW Push-to-talk		●	○	○
Transcribe		●	○	○
Multiflex Tinnitus Pro		●	○	○
ALL NEW Auracast™ assistant		●	●	●
ALL NEW Google Fast Pair™		●	●	●
Smart Assistant & 2-Way Audio		●	●	●
Hearing Health & Activity Tracking		●	●	●
Fall Alert & Balance Assessment		●	●	●
Reminders & Self Check		●	●	●
Tap Control		●	●	●
TeleHear full remote fit		●	●	●
REM Target Match		●	●	●
Demo Stream		●	●	●

Wireless Accessories

ALL NEW



StarLink Edge LE Adapter
Part No. P00018678



StarLink Edge Remote Microphone
Part No. P00004321



StarLink Edge Table Microphone
Part No. P00004326



StarLink Edge TV Streamer
Part No. P00004315




StarLink Remote Control 2.0
Part No. P00002347

AURACAST


Charging Solutions

All chargers include **Turbo Charge**


Serve as a charging station and carrying case for your rechargeable hearing aids.



StarLink Premium Charger 3.0
w/ onboard battery
• RIC RT



StarLink Premium Charger 2.0
w/ onboard battery
• mRIC R



StarLink Premium Custom Charger 2.0
w/ onboard battery
• ITC/ITER

StarLink Charger 3.0
• RIC RT

StarLink Charger 2.0
• mRIC R

StarLink Custom Charger 2.0
• ITC/ITER

Compatible Apps



My Audibel



Hear Share



Scan QR code to check phone compatibility



Download on the App Store



GET IT ON Google Play

Tech Tier/Channels/Bands
24 / 20* / 16*



Waterproof**
Rechargeable

For more product details, visit AudibelPro.com



Faceplate Colors

- Pink
- Light Brown
- Medium Brown
- Chestnut*
- Dark Brown
- Black

Shell Colors

(IIC only available in Black)

- Clear
- Pink
- Light Brown
- Blue/Red
- Black

Charging Solutions

All chargers include
Turbo Charge




StarLink Premium Charger 2.0
w/ onboard battery

- CIC R NW

StarLink Custom Charger 2.0

- CIC R NW

Fitting Ranges



Tech Tiers | Channels | Bands **24** **20*** **16***

Edge Mode+ (CIC styles only)

User Control ● ○ ○

Speech Optimization

Consonant Brightness ● Maximum ● Mild ● Mild

Sound Manager

Quiet		Maximum expansion	Moderate expansion	Mild expansion
Machine Noise Reduction		Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction
Speech in Noise Reduction		Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction
Transient Noise Reduction		Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 6 dB of reduction
Wind Reduction		Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction

Additional Key Features

Multiflex Tinnitus Pro with new tinnitus maskers ● ● ○
 REM Target Match ● ● ●
 Waterproof** ● ● ●



Based on the industry's most advanced processor and powered by our industry-leading Neuro Sound Technology, Audibel Signature Series deliver the size, performance and personalized fit that patients can only get from you.

*CIC styles only. **Up to 1m. Rechargeable style only.

Quality in... quality out

Your professional expertise is critical when selecting the best solution for a patient. Many factors and considerations come into play when determining if a patient is a custom hearing aid candidate, including:

- Audiogram
- Dexterity
- Ear anatomy and history
- Lifestyle

Ultimately, **the patient has come to you** to weigh all considerations—physical, audiological and lifestyle priorities—and recommend a solution that delivers both success and satisfaction.

A superior custom product fitting starts with a high quality ear impression. Making a perfect impression requires proper technique, attention to detail, and repetition.

As the custom hearing aid experts, we know a thing or two about making accurate, precision-tested impressions —and how good impressions can also be used to make earmolds, custom cases, hearing protection solutions and more.



Scan the QR code or go to AudibelPro.com for support materials, how-to videos, and more.

Developed with help from Audibel's most experienced impression takers, these tools provide best practice tips perfected from shooting thousands of impressions in our Center for Excellence.

VITALITY^{AI}

Tech Tier/Channels/Bands
24 / 20 / 16

AURACAST

CROS COMPATIBLE



mRIC R



RIC RT



RIC 312



- Waterproof*
- Rechargeable
- Wireless
- Telecoil

*Up to 1m.

Available Colors

- Standard Colors
- White
 - Tech Black
 - Graphite Gray
 - Silver
 - Chestnut
 - Caramel
 - Beige

SnapFit 2.0 Stock and Custom Cased Receivers



Comfort Bud Open
5 mm, 6 mm, 7 mm,
9 mm, 11 mm



Comfort Bud Occluded
4 mm, 5 mm, 6 mm,
7 mm, 9 mm, 11 mm



Power Dome
8-10 mm, 10-12 mm,
12-14 mm




Custom Cased Receiver
Custom shell with
embedded receiver



**Stock Silicone
Helix RIC Earmold**





**Stock Hollow
Silicone RIC Earmold**



**Stock Silicone
RIC Earmold**

Fitting Ranges



Visit AudibelPro.com/VitalityAI for full technical data sheets.



Powered by the G2 Neuro Processor with the neural processing unit (NPU) integrated into the chip, and Bluetooth® Low Energy (LE) Audio streaming protocol, Audibel's Vitality AI advanced processing is able to provide so much more than all-day battery life and waterproof protection – allowing patients to live easier, more connected lives.

For more Vitality AI details,
visit AudibelPro.com/VitalityAI



VITALITY^{AI}

Tech Tiers | Channels | Bands

24


20

16


Edge Mode+				
Automatic User Control Personalization App Control	●●●●	○●●●	○●●●	○●●●
Speech Optimization				
Consonant Brightness Voice AI*	●● Maximum	●○ Mild	●○ Mild	○● Mild
Sound Manager				
Quiet	●●●● Maximum expansion	●●●○ Moderate expansion	●●○○ Mild expansion	
Machine Noise Reduction	●●●● Up to 22 dB of reduction	●●●○ Up to 12 dB of reduction	●●○○ Up to 7 dB of reduction	
Speech in Noise Reduction	●●●● Up to 22 dB of reduction	●●●○ Up to 10 dB of reduction	●●○○ Up to 8 dB of reduction	
Transient Noise Reduction	●●●● Up to 15 dB of reduction	●●●○ Up to 9 dB of reduction	●●○○ Up to 6 dB of reduction	
Wind Reduction	●●●● Up to 35 dB of reduction	●●●○ Up to 19 dB of reduction	●●○○ Up to 10 dB of reduction	
Additional Key Features				
Translate	●	○	○	
Transcribe	●	●	○	
Multiflex Tinnitus Pro	●	●	○	
ALL NEW Auracast™ assistant	●	●	●	
Smart Assistant & Find My Hearing Aids	●	●	●	
2-Way Audio	●	●	●	
Hearing Health & Activity Tracking	●	●	●	
Fall Alert & Balance Assessment	●	●	●	
Reminders & Self Check	●	●	●	
Tap Control	●	●	●	
TeleHear full remote fit	●	●	●	
REM Target Match	●	●	●	
Demo Stream	●	●	●	
Waterproof up to 1m**	●	●	●	

Wireless Accessories

ALL NEW




StarLink Edge LE Adapter
Part No. P00018678




StarLink Edge Remote Microphone
Part No. P00004321



StarLink Edge Table Microphone
Part No. P00004326



StarLink Edge TV Streamer
Part No. P00004315




StarLink Remote Control 2.0
Part No. P00002347

AURACAST


Charging Solutions

All chargers include **Turbo Charge**




StarLink Premium Charger 3.0
w/ onboard battery

- RIC RT



StarLink Premium Charger 2.0
w/ onboard battery

- mRIC R



StarLink Premium Custom Charger 2.0
w/ onboard battery

- ITC/ITER

StarLink Charger 3.0


- RIC RT

StarLink Charger 2.0


- mRIC R

StarLink Custom Charger 2.0

- ITC/ITER



My Audibel



Hear Share



Scan QR code to check phone compatibility



Features noted above in green are available within My Audibel app.

*For iOS users only. **Rechargeable styles only.

INTRODUCING

A SERIES AI

Tech Tier/Channels/Bands
24* / 20* / 16* / 12 / 10



- Waterproof**
- Rechargeable
- Wireless
- Telecoil

*BTE styles only.
 **Up to 1m.

Available Colors


Standard Colors

- Graphite Gray
- Silver
- Chestnut
- Beige

SnapFit 2.0 Stock and Custom Cased Receivers



Fitting Ranges



Visit [AudibelPro.com](https://www.audibelpro.com) for full technical data sheets.



Powered by our industry-leading Neuro Sound Technology, A Series AI offers an easy, affordable, and immediate way to improve patient hearing and reconnect them with the world around them.







For more product details, visit [AudibelPro.com](https://www.audibelpro.com)



A SERIES^{AI}

Tech Tiers | Channels | Bands

24* **20*** **16*** **12** **10**

Speech Optimization					
Ear-to-Ear Technology	●	●	○	○	○
Frequency Lowering	●	●	●	●	●
Sound Manager					
Directionality 	Adaptive	Adaptive	Adaptive	Dynamic	Dynamic
Transient Noise 	●	○	○	○	○
Quiet 	Maximum expansion	●	●	●	●
Machine Noise Reduction 	Up to 22 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction
Speech in Noise Reduction 	Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction
Wind Reduction 	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction	Up to 10 dB of reduction	Up to 10 dB of reduction
Additional Key Features					
Edge Mode+	●	●	○	○	○
Reminders	●	●	○	○	○
Multiflex Tinnitus Technology	●	[white, pink, & ocean noise only]	[white, pink, & ocean noise only]	[white noise only]	○
Find My Hearing Aids	●	●	●	●	●
Self Check	●	●	●	●	●
TeleHear full remote fit	●	●	●	●	●
REM Target Match	●	●	●	●	●
Apple Watch Controls**	●	●	●	●	●

Wireless Accessories



StarLink Remote Control 2.0
Part No. P00002347



StarLink Mini Remote Microphone
Part No. 43121-000



StarLink Remote Microphone +
Part No. 43120-000



StarLink Edge Table Microphone
Part No. 43123-610



StarLink TV Streamer
Part No. 43119-000

Charging Solutions

All chargers include **Turbo Charge**



StarLink Premium Charger 3.0
w/ onboard battery
• RIC R



StarLink Charger 3.0
w/ onboard battery
• RIC R

Compatible Apps



My Audibel



Scan to check phone compatibility





Features noted above in green are available within My Audibel app.

*BTE styles only. **For iOS users only.

©2026 Starkey Laboratories, Inc. All Rights Reserved. P00008881 7/26 MISC5919-06-EN-AB

My Audibel App

Easy as it gets. Powerful as you need.



MyAudibel

Focus on benefits

Aris AI, Vitality AI, and A Series AI wireless hearing aids can connect via Bluetooth® technology to the My Audibel app.

- **Hear better anywhere**

Make instant adjustments for clearer speech and reduced noise—even in the toughest listening environments.

- **Get tips, tools, and resources**

Adapt to hearing aids faster with quick how-to's, tips, and videos. All tailored for where the hearing aid wearer is on their hearing journey.




- **Monitor health and wellness***

Monitor daily hearing and physical activity to help stay healthy. Set and achieve goals for hearing, steps, exercise, and beyond.


- **All from the palm of your hand or wrist**

Change programs and adjust volume quickly and easily on a compatible Apple Watch. Visit the Learn section within My Audibel for installation instructions.

How to download

1. Open the **App Store**  on Apple devices or **Google Play**  on Android devices.
2. In the search box, type  **My Audibel**.
3. Tap **GET** on Apple devices or **INSTALL** on Android devices. This will download the app to the compatible device.





A new dimension in connectivity

All-new **Auracast™ assistant** and **Google Fast Pair™** technology make connectivity so convenient, it's automatic.

ALL NEW

To learn more visit
AudibelPro.com



*Compatible with Aris AI and Vitality AI products only.

Other Apps


Enhance the capabilities of our hearing aids and simplify daily life with our lineup of popular hearing apps.



SoundCheck Live

Helps anyone gain awareness of noise exposure and the importance of protecting hearing by:

- Measuring sound levels around them using the built-in sound level meter.
- Learning about environmental noise's impact on hearing and how hearing protection devices like SoundGear can help.



Hear Share

Works with My Audibel app to allow wearers of Aris AI and Vitality AI to:

- Share information about hearing aid usage, social and physical activity progress, safety, and more—with the people they choose.
- Provide peace of mind for patients and their caregivers.



Hearing Protection

Be your community's full-service hearing care expert by offering hearing protection products from SoundGear. SoundGear products come in a variety of electronic and non-electronic solutions.



Electronic Solutions



PHANTOM Custom Fit Electronic

World's first custom, rechargeable, and Bluetooth® compatible hearing aid and protection device. Up to 30 dB gain.



SHIELD

Small, light and ready to take everywhere delivering dynamic digital hearing protection anytime you need it.

 **93 dB** | Part No. P00001867

 **80 dB** | Part No. P00001866



RACING RECEIVERS

Custom design allows for maximum comfort and durability for amateur and professional drivers to experience quality one-way listening communication.

	Top Seller	
	Phantom	Shield
Fit	Custom	Universal
Rechargeable	Yes	No
NRR	22 dB	26 dB
Warranty	1 year	1 year
Color Options	●	
Lanyard	●	●
Bluetooth	●	
Digital Modes	●	

Non-electronic Solutions



Filtered Earplug HS Silicone

Referred to as musician or industrial earplugs, these earplugs attenuate across all frequencies.



Solid Earplug HS Silicone

Outstanding noise protection in a variety of recreational and occupational environments.



High-Frequency Noise Filter HS Silicone

Referred to as "Hocks", these earplugs contain a patented acoustical chamber and filter that attenuate sound to appropriate levels.

Additional SoundGear Earmolds:

Newsreader/Teleprompter mold
Digital SLS/HS Silicone

Bluetooth mold
HS Silicone

Pilot headset with boom
Digital SLS

Sleeping earplug
HS Silicone

Anesthesiologist mold
Digital SLS

Swim earplug
Floatable Silicone

For a full selection of custom options, please contact Customer Service.

Personal Audio

For musicians, music aficionados, or professionals looking for a custom audio solution, Audibel has options to maximize their hearing experience.

Tunz Duo S | Tunz Duo | Tunz Solo

Additional Items

SoundGear Custom Charger for Phantom devices
Part No. 43059-610

Replacement Tri-flanges for Shield
Size 1: P00009811 (4pk)
Size 2: P00009812 (4pk)

Foam Buds Assorted for Shield
Part No. P00009810 (6pk)

Replacement filters
Filtered earplug: 10 dB, 17 dB, 26 dB, Impulse, Total block
High-Frequency (Hock): Part No. 50293-000

SoundGear Black Lanyard
Variety of lengths available
SoundGear Zipper Pouch
Part No. 91289-007

Hearing Aid Earmolds

As the leader in customization, Audibel understands offering your patients an expansive portfolio of hearing care solutions specific to their needs is key.



AUDIBEL

Material Options

Digital SLS Acrylic	<ul style="list-style-type: none"> • RIC • BTE & Thin Tube • Specialty* 	<ul style="list-style-type: none"> • One color 	Most Rigid 90 Shore
High Strength Silicone 60	<ul style="list-style-type: none"> • RIC • BTE & Thin Tube • Specialty* 	<ul style="list-style-type: none"> • Swirl (max. 3 colors) • Glitter option 	60 Shore
High Strength Silicone 40	<ul style="list-style-type: none"> • RIC • BTE & Thin Tube • Specialty* 	<ul style="list-style-type: none"> • Swirl (max. 3 colors) • Glitter option 	40 Shore
Floatable Silicone	<ul style="list-style-type: none"> • Swim earplug (SoundGear only) 	<ul style="list-style-type: none"> • Swirl (max. 3 colors) 	25 Shore Softest

Venting Guide

This guide is applicable for all hearing aid earmolds.** Vent size is limited by ear anatomy.

Vent Name	Approximate Size
Small	1.2 mm
Medium	1.8 mm
Large	2.55 mm
Select-A-Vent	3.37 mm
Extra Large	As large as canal allows

Specialty Solutions



HA Listening Test Mold
Digital SLS

A shell or skeleton style earmold with a tubing attachment for listening to ITE/ITC/CIC/RIC style and BTE hearing aids for listen checks.




Stent Mold
Digital SLS

A hollowed shell of canal style earmold to hold the ear canal open.

RIC Solutions

RIC Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Available for both SnapFit 2.0 and SnapFit Receivers



Custom Cased or Absolute Power (AP) Receivers Digital SLS


Cable Lengths:
1, 2, 3, 4, 5

Gain:
L, M, P, UP

Receivers:
SnapFit 2.0
(ten-pin connector)

SnapFit
(six-pin connector)

Original
(two-pin connector)



*Only available in clear

Pro tip: When measuring the receiver wire for a Custom Cased or AP receiver, report the number that correlates with the top of the entrance to the ear canal.

Color Options*

Our earmolds come in a variety of color options. Some colors appear differently in certain color families based on material type. **Swirl combinations are also available using up to three colors.**




Pro tip: Adding White as one of the three colors within a color swirl combination will help separate and make the other color choices "pop." Avoid combining light colors with "clear" or combining similar colors.

BTE Solutions

Thin Tube Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Tubing Lengths:


- 1
- 2
- 3, 3+
- 4, 4+
- 5, 5+



Traditional Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Tubing:

- 13 Regular
- 13 Heavy
- 13 Extra Heavy
- 13 Dry Tubing
- 13 Regular Snap Connect
- 13 Heavy Snap Connect



Business Solutions

In addition to our industry-leading products, we offer a full suite of business solutions designed to keep your practice thriving.



AUDIBEL

Business Equipment

State-of-the-art diagnostic equipment is essential to providing the best patient experience at your practice. With the help of our Business Equipment Services, it's never been easier to find and maintain the right solutions for all your testing and fitting needs— all in one place.

Equipment includes (but not limited to):

- **Audiometers**
- **Tympanometers**
- **Fitting equipment**
- **Combination units**
- **Video otoscopes**
- **Supplies**
- **Sound rooms/enclosures**

All Make Repair

Our All Make Repair department services hearing aids that are damaged and out of manufacturer's warranty. This even includes hearing aids by manufacturers no longer in business. **Our team is here to help with the following:**

- **All Make Repair estimates**
- **Reconditioned aids orders**
- **Miscellaneous parts and accessory orders**
- **Analog aids**
- **Body aids**
- **Bone conduction aids**
- **All Make Salvage Program**
- **Hearing Aid Recycle Program**

Batteries

Not only does every battery purchase deliver a positive jolt to your bottom line, each visit to your office for batteries is one more chance to build loyalty and generate referrals from your patients.

- **Drive retention**
Returning patients are more likely to purchase upgrades, warranties and accessories from you due to the ongoing care and service that you provide.
- **Build loyalty**
Superior aftercare will lock in your existing patients rather than allowing them to go to your competition for their future hearing healthcare needs.
- **Boost prospects**
Use battery promotions to attract new patients to your practice.
- **Protect patients' investment**
High quality batteries protect your patient's investment and enhance your patient/provider relationship.
- **Make the right selection**
Choose from long-lasting customized or branded batteries, depending on your business goals.

THE JOY OF HEARING

Audibel logo, Audibel, TeleHear, Audibel Vitality, Edge Mode, Multiflex Tinnitus, StarLink, Aris, DNN 360, TeleHear, HearShare, Sound [+ checkmark], Thrive, SoundGear logo, SoundGear, Phantom, A Series, Shield, The Joy of Hearing and Tunz are trademarks of Starkey Laboratories, Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by the Bluetooth SIG. Any use of such marks by Starkey is under license. Other trademarks and trade names are those of their respective owners.

Apple Watch is a trademark of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. "Made for iPod," "Made for iPhone" and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, iPad or Apple Watch may affect wireless performance. Google Play Store and Fast Pair are trademarks of Google, LLC.