

# Starkey's Complete Style Portfolio

Effortless selection of any style for every patient.



2.4 GHz HEARING AIDS

Feature/Product Style	ITC/HS R and ITE R	RIC R	BTE R	IIC NW*	CIC	RIC 312	micro RIC 312	BTE 13	Power Plus BTE 13
Battery Type	Rechargeable Lithium-ion	Rechargeable Lithium-ion	Rechargeable Lithium-ion	10	312	312	312	13	13
Telecoil	• (ITE R only)	•	•					•	•
CROS System		•				•		•	
Tinnitus Technology	•	•	•	•	•	•	•	•	•

\*IIC NW excluded from Thrive 2.4 GHz platform

2.4 GHz ACCESSORIES



## Table Microphone

With our 2.4 GHz Table Microphone, users can enjoy group settings like family gatherings or noisy restaurants. Eight built-in microphones detect the primary speaker's voice and stream it directly to their hearing aids.



## TV Streamer

Stream audio from the television or other electronic audio source directly to their hearing aids. It offers excellent sound quality, is easy to use and supports both analog and digital input sources.



## Remote Mic Plus

Use as a remote microphone, or easily stream binaurally from four different audio sources: Bluetooth®, Telecoil, FM, or Line-In.



## Mini Remote Microphone

Enjoy one-on-one conversations in noisy environments with our small, easy-to-use Mini Remote Microphone. Clip it onto the clothing of the person being spoken to, or use it as a TV streamer by placing it near the speaker.



## Remote Control

Our Remote Control includes updated features to control memory and volume, mute hearing aids and turn other special features on and off.



## Premium Standard Charger & Custom Charger

Both chargers hold enough power to provide portable charging three times. LED indicators let patients know when their hearing aids are fully charged, and hearing aids turn on automatically when removed from the charger.



## Mini Turbo Charger

The Mini Turbo Charger is a pocket-size lithium-ion charging option providing a 3.5 hour charge in just 7 minutes or a full charge in under 3 hours. This portable charger will recharge a pair of RIC R or BTE R hearing aids up to four times without a power connection.



# Thrive Hearing Control App



## Hearing Performance & Health Monitoring Technology Features:

Available with **Evolv AI** products



### Edge Mode

With just a double tap, Edge Mode instantly conducts an AI-based analysis of the listening environment and makes smart, immediate adjustments.



### Engagement and Activity Tracking

Monitors vigorous physical activity and Health monitoring features like active listening.



### Fall Alert

With integrated sensors, Fall Alert can detect falls and send alert messages to selected contacts.



### Find My Phone

Using hearing aids and Thrive Assistant\*\*, you can now locate your lost or misplaced phone by having it ring. No other hearing aid offers this capability.



### Mask Mode

A custom memory that boosts speech audibility.



### Self Check

Allows you to perform a diagnostic test of your hearing aid components.



### Thrive Assistant

You can conveniently troubleshoot hearing aid and accessory issues and ask a question like, "How do I customize my hearing aid settings?"



### Personalized Control

Make personalized, real-time adjustments.



### 2 Way Audio

Now, you can take and end calls directly from your hearing aids instead of your Apple® device. Plus, hearing aid microphones pick up and stream your voice directly back to iPhone and iPad\*, allowing you to seamlessly talk hands-free.



### Language Translation

Easily communicate with people who speak other languages via the easy-to-use translation tool in your Thrive app. (Requires an internet connection.)



### Tap Control

Tap your ear twice to start & stop streaming audio from a 2.4 GHz accessory.



### TeleHear®

Remote Programming allows for virtual adjustments to hearing aids.



### Thrive Assistant

Use voice commands\*\* to get answers to questions, set reminders and control hearing aids.



### Remote Control

Change volume and switch memories with your smartphone.



### Personalized Memories

Create custom memories, and geotag memories so adjustments occur automatically when you return to that location.



### Auto Experience Manager

Automatically adjusts loudness over time to help you transition your new hearing experience comfortably.



### Phone Calls

Answer calls and hear the conversation streamed directly to your hearing aids using compatible devices.



### Transcribe

Conversations are transcribed so you can read what's being said into text that can be saved or shared.



### Voice AI\*

This always-learning speech enhancement feature uses an iOS smartphone to deliver an unprecedented experience for those with moderate loss of 50 dB HL or greater.



### Find My Hearing Aids

Easily locate lost hearing aids.



### Adaptive Car Mode

Hearing aids automatically change to reduce the annoying sounds of driving.



### Reminders

Personalized or on-screen reminders for daily tasks such as taking medications.



### Siri Integration

Have Siri® read texts and emails directly to your hearing aids (available with Apple devices).

\* Requires an iOS device. \*\*Voice-activated commands and reminders are exclusive to Evolv AI.

The app may look slightly different depending on your phone.

Starkey, Starkey logo, Thrive, Thrive logo, Edge Mode and Evolv are registered trademarks of Starkey Laboratories, Inc. Evolv logo is a trademark of Starkey Laboratories, Inc. Google Play is a trademark of Google, LLC. | The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks is under license. | Siri and Apple are registered trademarks of Apple Inc.

