


Product Portfolio Guide



Future Forward

By continuously pushing the edge of what's possible, Audibel will always be 10 years ahead—with groundbreaking innovations, award-winning products and services, and a **complete ecosystem of solutions** meticulously integrated to help you effortlessly provide joyful hearing for all.





Available Colors




Helix RIC Earmold

SnapFit 2.0 Stock and Custom Cased Receivers

Silicone RIC Earmold

Fitting Ranges



Visit AudibelPro.com/VitalityAI for full technical data sheets.



Powered by the G2 Neuro Processor with the industry's only neural processing unit (NPU) integrated into the chip, and an all-new Bluetooth® Low Energy (LE) Audio streaming protocol, Audibel Vitality AI's advanced processing is able to provide so much more than all-day battery life and waterproof protection – allowing patients to live easier, more connected lives.

For more Vitality AI details, visit AudibelPro.com/VitalityAI

RIC Earmold



VITALITYA



Wireless Accessories





StarLink Edge Table Microphone Edge TV Streamer Part No. P00004326



StarLink **Remote Control 2.0** Part No. P00002347



w/onboard battery



StarLink Premium Custom Charger 2.0

• ITC/ITE R: Part No. P00009158

StarLink Custom Charger 2.0***

• ITC/ITE R: Part No. P00009157





StarLink Premium Charger 2.0 w/ onboard battery

• mRIC R: Part No. P00009110 • RIC RT: Part No. P00009144

StarLink Charger 2.0***


• mRIC R: Part No. P00009108 • RIC RT: Part No. P00009109

***Default selection with rechargeable hearing aid order.


Compatible Apps



Features noted above in green are available within My Audibel app.



Neuro Sound Technology, Audibel Signature Series deliver the size, performance and personalized fit that patients can only get from you.



Quality in... quality out

Your professional expertise is critical when selecting the best solution for a patient. Many factors and considerations come into play when determining if a patient is a custom hearing aid candidate, including:

- Audiogram
- Dexterity
- Ear anatomy and history
- Lifestyle

Ultimately, **the patient has come to you** to weigh all considerations—physical, audiological and patient priorities—and recommend a solution that delivers both success and satisfaction.

A superior custom product fitting starts with a high quality ear impression. Making a perfect impression requires proper technique, attention to detail, and repetition.

As the custom hearing aid experts, we know a thing or two about making accurate, precision-tested impressions — and how good impressions can also be used to make earmolds, custom cases, hearing protection solutions and more.



Scan the QR code or go to AudibelPro.com for support materials, how-to videos, and more.

Developed with help from Audibel's most experienced impression takers, these tools share best practice tips perfected from shooting thousands of impressions in our Center for Excellence.




Available Colors

SnapFit 2.0 Stock and Custom Cased Receivers



Fitting Ranges



Visit AudibelPro.com for full technical data sheets.



Intrigue AI marked a new beginning for hearing technology by introducing the industry to the most advanced processor technology with a unique onboard Deep Neural Network (DNN) accelerator. This Neuro Sound Technology in Intrigue AI is designed to deliver sound so transparent and true to life, patients will forget they're wearing hearing aids.

For more product details, visit AudibelPro.com



INTRICUEA

Tech Tiers Channels Bands	24	20	16	12
	24	20	10	12
Edge Mode+ Automatic User Control Personalization App Control	• • •		0 0 0	0 0 0
Speech Optimization				
Consonant Brightness Voice AI*	Maximum	Mild	Mild	0
Sound Manager				
Quiet 🛞	Maximum expansion	Moderate expansion	Mild expansion	• O O Mild expansion
Machine Noise Reduction $6 = 6$	Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction	• O O O Up to 7 dB of reduction
Speech in Noise Reduction	Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction	• • • • • Up to 8 dB of reduction
Transient Noise Reduction ())	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 6 dB of reduction	• • • • • Up to 6 dB of reduction
Wind Reduction	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction	• • • • • Up to 10 dB of reduction
Additional Key Features				-
Translate Transcribe Multiflex Tinnitus Pro Smart Assistant Find My Hearing Aids 2-Way Audio Hearing Health & Activity Tracking Fall Alert Reminders & Self Check Auto Sleep & Tap Control TeleHear full remote fit REM Target Match Waterproof up to 1m** Apple Watch Controls*				

Wireless Accessories







StarLink Premium Custom Charger 2.0

• ITC/ITE R: Part No. P00009158

StarLink Custom Charger 2.0***

• ITC/ITE R: Part No. P00009157

w/onboard battery

StarLink



StarLink StarLink Mini Remote Control 2.0 Remote Microphone Remote Microphone + Table Microphone Part No. P00002347 Part No. 43121-000 Part No. 43120-000 Part No. 43123-610

StarLink



Charging Solutions





StarLink Premium Charger 2.0 w/ onboard battery

• mRIC R: Part No. P00009110 • RIC RT: Part No. P00009144

StarLink Charger 2.0***

• mRIC R: Part No. P00009108 • RIC RT: Part No. P00009109

***Default selection with rechargeable hearing aid order.


Compatible Apps



Features noted above in green are available within My Audibel app.



Available Colors



Stock Solutions

Comfort Bud Open

5 mm, 6 mm, 7 mm,

9 mm, 11 mm






Power Dome

12-14 mm



Stock Thin Tube Solution 8-10 mm, 10-12 mm.

Fitting Ranges




Power Plus BTE 13

Comfort Bud Occluded

4 mm, 5 mm, 6 mm,

7 mm, 9 mm, 11 mm



Visit AudibelPro.com for full technical data sheets.

Arc Al offers a complete line of BTE (Behind-The-Ear) styles loaded with innovative technology and clear sound that seamlessly connect your patients to the life they love.

For more product details, visit AudibelPro.com



Tech Tiers	2400	2000	1600	1200	1000
Channels Bands	24	20	16	12	10
Edge Mode					
User Control App Control	•			0	0
Speech Optimization					
Spatial Speech Enhancement Ear-to-Ear Technology Voice Al*	•			0000	000
Sound Manager		÷		-	
Auto Music	Maximum expansion	Moderate expansion	Mild expansion	Mild expansion	Mild expansion
Quiet 🔗	Maximum expansion	Moderate expansion	Mild expansion	Mild expansion	Mild expansion
Machine Noise Reduction	Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction
Speech in Noise Reduction	Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB	Up to 8 dB of reduction	Up to 8 dB of reduction
Transient Noise Reduction))	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction
Wind Reduction	Up to 32 dB of reduction	Up to 17 dB of reduction	Up to 9 dB of reduction	Up to 9 dB of reduction	Up to 9 dB of reduction
Additional Key Features					
Translate Multiflex Tinnitus Pro Transcribe Thrive Assistant & Find My Hearing Aids CROS System** Engagement Tracking	•				0 0 0 0 0
Activity Tracking Fall Alert 2-Way Audio Reminders & Self Check	•	•			
Auto On/Off & Tap Control TeleHear full remote fit REM Target Match	•	•	•	• •	•

Wireless Accessories







StarLink StarLink Mini Remote Control 2.0 Remote Microphone Remote Microphone + Part No. P00002347 Part No. 43121-000

StarLink Part No. 43120-000





StarLink StarLink **Table Microphone TV Streamer** Part No. 43123-610 Part No. 43119-000

Charging Solutions



Standard Premium Charger w/onboard battery • BTE R: Part No. P00010092

Standard Charger 2.0*** • BTE R: Part No. P00010091

***Default selection with rechargeable hearing aid order.

Compatible Apps





Scan to check phone compatibility



GET IT ON **Google** Play

Features noted above in green are available within Thrive app. *For iOS users only. **BTE 13 style only.

Hearing Protection

Be your community's full-service hearing care expert by offering hearing protection products from SoundGear. SoundGear products come in a variety of electronic and non-electronic solutions.



Electronic Solutions



PHANTOM **Custom Fit Electronic**

World's first custom, rechargeable, and Bluetooth[®] compatible hearing aid and protection device. Up to 30 dB gain.



ALL-NEW SHIELD

Small, light and ready to take everywhere delivering dynamic digital hearing protection anytime you need it.



93 dB | Part No. P00001867

80 dB | Part No. P00001866

Non-electronic Solutions



Filtered Earplug HS Silicone

Referred to as musician or industrial earplugs, these earplugs attenuate across all frequencies.



Solid Earplug HS Silicone

Outstanding noise protection in a variety of recreational and occupational environments.



High-Frequency Noise Filter HS Silicone

Referred to as "Hocks", these earplugs contain a patented acoustical chamber and filter that attenuate sound to appropriate levels.

Silver **Custom Fit Electronic**

ITE custom fit electronic molded to each individual's unique ear. Up to 30 dB gain.

Racing Receiver

Custom design allows for maximum comfort and durability for amateur and professional drivers to experience quality one-way listening communication.

	Top Seller		
	Phantom	Silver	Shield
Fit	Custom	Custom	Universal
Rechargeable	Yes	No	No
NRR	22 dB	24 dB	26 dB
Warranty	1 year	1 year	1 year
Color Options	•	٠	
Lanyard	•	٠	٠
Bluetooth	•		
Digital Modes	•		

Additional SoundGear earmolds:

Newscaster/Teleprompter Mold	Bluetooth mold
Digital SLS/HS Silicone	HS Silicone

Pilot Headset with Boom Digital SLS

Sleeping earplug HS Silicone

Anesthesiologist mold Digital SLS

Swim earplug Floatable Silicone

For a full selection of custom options. please contact Customer Service.

Personal Audio

For musicians, music aficionados, or professionals looking for a custom audio solution, Starkey has options to maximize their hearing experience.

Tunz Duo S | Tunz Duo | Tunz Solo

Additional Items

SoundGear Custom Charger for Phantom devices Part No 43059-610

Replacement Tri-Flanges for Shield Size 1: P00009811 (4pk) Size 2: P00009812 (4pk)

Foam Buds Assorted for Shield

Part No. P00009810 (6pk)

Replacement filters Filtered earplug: 10 dB, 17 dB, 26 dB, Impulse, Total block **High-Frequency** (Hock): Part No. 50293-000

SoundGear Black Lanyard Variety of lengths available

SoundGear Zipper Pouch Part No 91289-007

Hearing Aid Earmolds

As the leader in customization, Audibel understands offering your patients an expansive portfolio of hearing care solutions specific to their needs is key.



Canal

Open

Canal

Lock

Material Options

Digital SLS	• RIC • BTF & Thin Tube	• One color	Most Rigid			
Acrylic	· Specialty'		90 Shore			
High Strength Silicone 60	 RIC BTE & Thin Tube Specialty* 	 Swirl (max. 3 colors) Glitter option 	60 Shore			
High Strength Silicone 40	 RIC BTE & Thin Tube Specialty* 	 Swirl (max. 3 colors) Glitter option 	40 Shore			
Floatable Silicone	 Swim earplug (SoundGear only) 	• Swirl (max. 3 colors)	25 Shore Softest			
			Juliest			

Venting Guide

This guide is applicable for all hearing aid earmolds. Vent size is limited by ear anatomy.

Vent Name	Approximate Size
Pressure Vent	0.8 mm
Small	1.2 mm
Medium	1.8 mm
Large	2.55 mm
Select-A-Vent	3.37 mm
Extra Large	As large as canal allows

Specialty Solutions



HA Listening Test Mold Digital SLS

Stent Mold

A shell or skeleton style earmold with a tubing attachment for listening to ITE/ITC/CIC/RIC style and BTE hearing aids for listen checks.



Digital SLS A hollowed shell of canal style earmold to hold the ear canal open.

RIC Solutions

RIC Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40 Thin Tube Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40 Available for both SnapFit 2.0 and SnapFit Receivers **Tubing Lengths:** 1 2 3,3+ 4,4+ RIC RIC Half Shell Canal RIC RIC RIC CIC Shell 5,5+ Hollow Canal Canal Solid Canal Lock Skeleton Helix RIC Canal Lock Canal Skeleton Helix Helix Occluded Lock Open Custom Cased or Absolute Power (AP) Receivers Digital SLS Traditional Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40 Cable Lengths: *Only available in clear Tubing: 1, 2, 3, 4, 5 13 Regular Pro tip: When Gain: measuring the 13 Heavy L, M, P, UP receiver wire for 13 Extra Heavy a Custom Cased CIC Shell Half Canal RIC RIC or AP receiver. Receivers: Shell Custom Cased / Custom Cased / 13 Dry Tubing report the number SnapFit 2.0 AP Canal AP Canal Lock that correlates 13 Regular Snap (ten-pin connector) with the top of the Connect entrance to the SnapFit ear canal. 13 Heavy Snap (six-pin connector) Connect Skeleton Original (two-pin connector)

BTE Solutions

Color Options^{*}

Our earmolds come in a variety of color options. Some colors appear differently in certain color families based on material type. **Swirl combinations are also available using up to three colors.**

Digital S	SLS																	
												Pro Tip: /						
Clear	Rose Transp.	Tan	Brown	Red	Blue	Green	Yellow	Purple	Orange	Black		color swi other col			ill help s	eparate	and ma	ike the
HS Silic	one 60		1	1								Color cor					, ,	
Clear	Rose	Brown	Blue	Green	Yellow	Orange	Black	Camou-				Light' o' combin	colors co ed with '		n Swirt d	or Light	colors	
oteal	Transp.	DIOWII	Diue	oreen	Tettow	orange	DIACK	flage**				 Similar 	colors c	ombined	l in Swirl			
HS Silic	one 40			1										_				
Clear	Rose Transp.	Light Brown	Brown	Red	Blue	Green	Yellow	White	Purple	Neon Pink	Orang	je Light Blue	Light Pink	Light Purple	Cham- pagne	Onys	Slate	Sterling
Glitter																		
Blue	Green	Pink	Red	Orange	Gold	Silver	Purple	Multi- Colored										
Floatab	le**																	
Rose Transp.	Tan	Red	Blue	Green	Yellow	White	Purple	Neon Pink	Orange	Light Blue	Blac	¢						

*Also applies to SoundGear non-electronic solutions. **Available for SoundGear only.

Business Solutions

In addition to our industry-leading products, we offer a full suite of business solutions designed to keep your practice thriving.



Business Equipment

State-of-the-art diagnostic equipment is essential to providing the best patient experience at your practice. With the help of our Business Equipment Services, it's never been easier to find and maintain the right solutions for all your testing and fitting needs— all in one place. Equipment includes (but not limited to):

- Audiometers
- Tympanometers
- Fitting equipment
- Combination units
- Video otoscopes
- Supplies
- Sound rooms/enclosures

All Make Repair

Our All Make Repair department services hearing aids that are damaged and out of manufacturer's warranty. This even includes hearing aids by manufacturers no longer in business. **Our team is here** to help with the following:

- All Make Repair estimates
- Reconditioned aids orders
- Miscellaneous parts and accessory orders
- Analog aids
- Body aids
- Bone conduction aids
- All Make Salvage Program
- Hearing Aid Recycle Program

Batteries

Not only does every battery purchase deliver a positive jolt to your bottom line, each visit to your office for batteries is one more chance to build loyalty and generate referrals from your patients.

Drive retention

Returning patients are more likely to purchase upgrades, warranties and accessories from you due to the ongoing care and service that you provide.

Build loyalty

Superior aftercare will lock in your existing patients rather than allowing them to go to your competition for their future hearing healthcare needs.

Boost prospects

Use battery promotions to attract new patients to your practice.

- **Protect patients' investment** High quality batteries protect your patient's investment and enhance your patient/provider relationship.
- Make the right selection

Choose from long-lasting customized or branded batteries, depending on your business goals.



Scan to purchase.

THE JOY OF HEARING

Audibel logo, Audibel, TeleHear, Pro Fit, Audibel Vitality, Edge Mode, Multiflex Tinnitus, StarLink, Audibel, Neuro Sound Technology, Hear Share, Intrigue, Thrive, SoundGear logo, SoundGear, Phantom, Arc, Shield, The Joy of Hearing and Tunz are trademarks of Starkey Laboratories, Inc.

The Bluetooth[®] word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast[™] word mark and logos are trademarks owned by the Bluetooth SIG. Any use of such marks by Starkey is under license. Other trademarks and trade names are those of their respective owners.

Apple Watch is a trademark of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. "Made for iPod," "Made for iPhone" and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, iPad or Apple Watch may affect wireless performance. Google Play Store is a trademark of Google, LLC.