


# **Product Portfolio Guide**

Fall 2024



# Future Forward

By continuously pushing the edge of what's possible, Audibel will always be 10 years ahead—with groundbreaking innovations, award-winning products and services, and a **complete ecosystem of solutions** meticulously integrated to help you effortlessly provide joyful hearing for all.











**♦** ⋒ **\*** 




Waterproof\* **Wireless** 

Rechargeable

\*Up to 1m.

#### **Available Colors**





**SnapFit 2.0 Stock and Custom Cased Receivers** 



Stock Silicone

Helix RIC Earmold

Comfort Bud Open Comfort Bud Occluded 5 mm, 6 mm, 7 mm, 4 mm, 5 mm, 6 mm, 9 mm, 11 mm





Stock Hollow Silicone RIC Earmold

**Power Dome** 8-10 mm, 10-12 mm. 12-14 mm





Custom Cased Receiver Custom shell with embedded receiver




Stock Silicone RIC Earmold

# **Fitting Ranges**








Visit AudibelPro.com/VitalityAI for full technical data sheets.



Powered by the G2 Neuro Processor with the industry's only neural processing unit (NPU) integrated into the chip, and an all-new Bluetooth® Low Energy (LE) Audio streaming protocol, Audibel Vitality Al's advanced processing is able to provide so much more than all-day battery life and waterproof protection - allowing patients to live easier, more connected lives.







#### **Wireless Accessories**







Remote Control 2.0 Part No. P00002347

**AURACAST** 

#### **Charging Solutions**

All chargers include

Turbo Charge —



StarLink Premium Charger 2.0 w/onboard battery

- mRIC R
- RIC RT



StarLink Premium Custom Charger 2.0 w/onboard battery

• ITC/ITE R

StarLink Charger 2.0\*\*\*

- mRIC R
- RIC RT

StarLink Custom Charger 2.0\*\*\*

• ITC/ITE R

# **Compatible Apps**







Scan to check phone compatibility





**Hear Share** 

<sup>\*\*\*</sup>Default selection with rechargeable hearing aid order.



# Nothing is smaller. Nothing sounds better.



Tech Tier/Channels/Bands **24 / 20\* / 16\*** 







6 4



For more product details, visit AudibelPro.com



#### **Faceplate Colors**

- Pink
- Light Brown
- Medium Brown
- Chestnut\*Dark Brown
- Black

#### **Shell Colors**

(IIC only available in Black)

- Clear
- Pink
- Light Brown
- Blue/Red
- Black

#### **Charging Solutions**

All chargers include

Turbo Charge




StarLink Premium Charger 2.0 w/onboard battery

• CIC RW

StarLink Custom Charger 2.0\*\*\*

• CIC R NW

# **Fitting Ranges**





Based on the industry's most advanced processor and powered by our industry-leading Neuro Sound Technology, Audibel Signature Series deliver the size, performance and personalized fit that patients can only get from you.

#### **Premium Advanced** Select 24 20 16\* Tech Tiers | Channels | Bands Edge Mode+ (CIC styles only) User Control Speech Optimization Consonant Brightness Mild Mild Maximum Sound Manager 00 Quiet

Speech in Noise Reduction

Machine Noise Reduction





Maximum

expansion

Up to 22 dB



Moderate

expansion

Up to 12 dB

of reduction

Up to 8 dB of reduction

Mild

expansion

000

Up to 7 dB

of reduction

Transient Noise Reduction



Up to 15 dB of reduction

Up to 9 dB of reduction

Up to 6 dB of reduction

Wind Reduction



Up to 35 dB of reduction

Up to 19 dB of reduction

Up to 10 dB of reduction

#### Additional Key Features

Multiflex Tinnitus Pro with new tinnitus maskers REM Target Match Waterproof\*\*





•



# Quality in... quality out

Your professional expertise is critical when selecting the best solution for a patient. Many factors and considerations come into play when determining if a patient is a custom hearing aid candidate, including:

- Audiogram
- Dexterity
- Ear anatomy and history
- Lifestyle

Ultimately, **the patient has come to you** to weigh all considerations—physical, audiological and patient priorities—and recommend a solution that delivers both success and satisfaction.

A superior custom product fitting starts with a high quality ear impression. Making a perfect impression requires proper technique, an attention to detail, and repetition.

As the custom hearing aid experts, we know a thing or two about making accurate, precision-tested impressions —and how good impressions can also be used to make earmolds, custom cases, hearing protection solutions and more.



Scan the QR code or go to AudibelPro.com for support materials, how-to videos, and more.

**Developed with help from Audibel's most experienced impression takers**, these tools share best practice tips perfected from shooting thousands of impressions in our Center for Excellence.

# INTRIGUEAL

Tech Tier/Channels/Bands 24 / 20 / 16 / 12





CIC

3



ITC R

**♦** ⋒ **\*** 



ITE R


**6 a 7** 



- Waterproof\*
- **Wireless**
- Rechargeable

\*Up to 1m.

#### **Available Colors**





Black

**SnapFit 2.0 Stock and Custom Cased Receivers** 









Comfort Bud Open 5 mm, 6 mm, 7 mm, 9 mm, 11 mm

4 mm, 5 mm, 6 mm, 7 mm, 9 mm, 11 mm

**Power Dome** 8-10 mm, 10-12 mm. 12-14 mm

Custom Cased Receiver Custom shell with embedded receiver



Stock Silicone Helix RIC Earmold





Stock Hollow Silicone RIC Earmold




Stock Silicone RIC Earmold

# **Fitting Ranges**







Visit AudibelPro.com for full technical data sheets.



Intrique AI marked a new beginning for hearing technology by introducing the industry to the most advanced processor technology with a unique onboard Deep Neural Network (DNN) accelerator. This Neuro Sound Technology in Intrigue AI is designed to deliver sound so transparent and true to life, patients will forget they're wearing hearing aids.



#### INTRIGUEAL Premium Advanced Select Basic 24 20 16 12 Tech Tiers | Channels | Bands Edge Mode+ Automatic User Control Personalization App Control Speech Optimization **Consonant Brightness** Mild Mild Maximum Voice AI\* Sound Manager 000 Quiet • 0 0 Maximum Moderate Mild Mild expansion expansion expansion expansion 000 • 0 0 0 Machine Noise Reduction Up to 22 dB Up to 12 dB Up to 7 dB Up to 7 dB of reduction of reduction of reduction of reduction • 0 0 0 Speech in Noise Reduction Up to 22 dB Up to 10 dB Up to 8 dB Up to 8 dB of reduction of reduction of reduction of reduction | -||›) **Transient Noise Reduction** • 0 0 0 Up to 15 dB Up to 9 dB Up to 6 dB Up to 6 dB of reduction of reduction of reduction of reduction • 0 0 0 Wind Reduction Up to 35 dB Up to 19 dB Up to 10 dB Up to 10 dB of reduction of reduction of reduction of reduction Additional Key Features Translate Transcribe Multiflex Tinnitus Pro Smart Assistant & Find My Hearing Aids 2-Way Audio Hearing Health & Activity Tracking Fall Alert

#### **Wireless Accessories**





StarLink Mini

Part No. P00002347

Remote Control 2.0 Remote Microphone Remote Microphone + Table Microphone Part No. 43121-000

Part No. 43120-000 Part No. 43123-610




StarLink TV Streamer Part No. 43119-000

# **Charging Solutions**

All chargers include Turbo Charge —





- mRIC R
- RIC RT

StarLink Charger 2.0\*\*\*

- mRIC R
- RIC RT



StarLink Premium Custom Charger 2.0 w/onboard battery

• ITC/ITE R

StarLink Custom Charger 2.0\*\*\*

• ITC/ITE R

# **Compatible Apps**







Scan to check phone compatibility





My Audibel

**Hear Share** 

Reminders & Self Check

TeleHear full remote fit

Waterproof up to 1m\*\*

Apple Watch Controls\*

**REM Target Match** 

Auto Sleep

Tap Control

<sup>\*\*\*</sup>Default selection with rechargeable hearing aid order.



Tech Tiers 2400 / 2000 / 1600 / 1200 / 1000





⋒ Wireless**∮** Rechargeable

#### **Available Colors**

Champagne

Slate

Sterling

Bronze

Black

Espresso

Bright White with Sterling

#### **Stock Solutions**



Comfort Bud Open 5 mm, 6 mm, 7 mm, 9 mm, 11 mm



Comfort Bud Occluded 4 mm, 5 mm, 6 mm, 7 mm, 9 mm, 11 mm




**Power Dome** 8-10 mm, 10-12 mm, 12-14 mm



Stock Thin Tube Solution

# **Fitting Ranges**





Arc <sup>AI</sup> Tech Tiers  Channels   Bands	Premium <b>2400</b> 24	Advanced 2000	Select <b>1600</b> 16	Basic <b>1200</b> 12	Entry 1000 10
Edge Mode					
User Control App Control	•	0	•	0	0
Speech Optimization					
Spatial Speech Enhancement Ear-to-Ear Technology Voice AI*	•	•	•	0 0	0 0
Sound Manager					
Auto Music	Maximum expansion	Moderate expansion	Mild expansion	Mild expansion	Mild expansion
Quiet 🚱	Maximum expansion	Moderate expansion	Mild expansion	Mild expansion	Mild expansion
Machine Noise Reduction	Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction
Speech in Noise Reduction	Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction
Transient Noise Reduction (1)	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction
Wind Reduction	Up to 32 dB of reduction	Up to 17 dB of reduction	Up to 9 dB of reduction	Up to 9 dB of reduction	Up to 9 dB of reduction
Additional Key Features					
Translate  Multiflex Tinnitus Pro Transcribe Thrive Assistant & Find My Hearing Aids CROS System** Engagement & Activity Tracking Fall Alert 2-Way Audio Reminders & Self Check Auto On/Off Tap Control TeleHear full remote fit REM Target Match					

#### **Wireless Accessories**







StarLink Remote Control 2.0 Remote Microphone Remote Microphone + Part No. P00002347 Part No. 43121-000 Part No. 43120-000







StarLink TV Streamer

# **Charging Solutions**



Standard Premium Charger w/onboard battery

• BTE R

Standard Charger 2.0\*\*\*

• BTE R

**HEAR** Share

Hear Share

# **Compatible Apps**







phone compatibility





<sup>\*\*\*</sup>Default selection with rechargeable hearing aid order.

# **Hearing Protection**

Be your community's full-service hearing care expert by offering hearing protection products from SoundGear. SoundGear products come in a variety of options for outdoor, work, or race.

FIND ME: eStore > Shop > Hearing Protection



#### **Electronic custom solutions**











# PHANTOM

#### **Custom Fit Electronic**

World's first custom, rechargeable, and Bluetooth® compatible hearing aid and protection device.



#### Silver **Custom Fit Electronic**

ITE custom fit electronic molded to each individual's unique ear.





#### **Racing Receiver**

Custom design allows for maximum comfort and durability for amateur and professional drivers to experience quality one-way listening communication.



#### **Instant Fit Electronic**

Uses: Hunting, Industry, Tactical Shooting, USPSA, 3-Gun, Sporting Clays, Trap and Skeet



Industrial • 8 dB gain Part No. 65720-001



Recreational • 15 dB gain Part No. 65720-000

#### Non-electronic custom solutions



#### **Mini Filtered Plugs**

Referred to as musician or industrial earplugs, these earplugs attenuate across all frequencies.



#### Solid Plugs

Outstanding noise protection in a variety of recreational and occupational environments.



#### **High-Frequency Noise Filters**

Referred to as "Hocks", these earplugs contain a patented acoustical chamber and filter that attenuate sound to appropriate levels in the workplace or on the range.

	Mini Filter Plugs		High-Frequency Noise Filters	
Fit	Custom	Custom	Custom	
NRR	Variable	26 dB	19-22 dB	
Shell Material	High Strength Silicone	High Strength Silicone	High Strength Silicone	
Color Options	•	•	•	
Lanyard	•	•	•	
Warranty	90 days	90 days	90 days	

# **Additional items**

SoundGear Custom Charger for Phantom devices Part No. 43059-610

Replacement sleeves for Instant Fit Electronics

Small (6pk): Part No. 52627-000 Large (6pk): Part No. 52628-000 Replacement Filters for Filtered Earplugs

10 dB, 17 dB, 26 dB, Impulse, Total block

SoundGear Black Lanyard Part No. 31274-000

SoundGear Zipper Pouch Part No. 91289-007

aficionados, or professionals looking for a custom audio solution, we have options to maximize their hearing experience.

Personal Audio: For musicians, music

Tunz Solo | Tunz Duo | Tunz Duo S

# **Earmolds**

As the leader in customization, Audibel understands offering your patients an expansive portfolio of hearing care solutions specific to their needs is key.

FIND ME: eStore > Shop > Earmolds



# **Material Options**

Digital SLS	· RIC	· One color	Most Rigid
Acrylic	· Specialty	· One cotor	90 Shore
High Strength Silicone 60	· RIC · BTE & Thin Tube · Specialty	Swirl 3 colors     Glitter option	60 Shore
High Strength Silicone 40	· RIC · BTE & Thin Tube · Specialty	Swirl 3 colors     Glitter option	40 Shore
Floatable Silicone	· RIC · BTE & Thin Tube · Specialty	· Swirl 3 colors	25 Shore Softest

# **Venting Guide**

This guide is applicable for all earmolds and earplugs. Vent size is limited by ear anatomy.

Vent Name	Approximate Size
Pressure Vent	0.8 mm
Small	1.2 mm
Medium	1.8 mm
Large	2.55 mm
Select-A-Vent	3.37 mm
Extra Large	As large as canal allows

# **Color Options**

Our earmolds come in a variety of color options. The latest color options can be found on the eStore. Some colors appear differently in certain color families based on material type. Swirl combinations are also available using up to three colors.

Pro Tip: Adding White as one of the three colors within a color swirl combination will help separate and make the other color choices "pop."

#### Color combinations to avoid:

- 'Light' colors combined in Swirl or 'Light' colors combined with 'Clear'
- · Similar colors combined in Swirl

#### **RIC Solutions**

#### RIC Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

#### Available for both SnapFit 2.0 and SnapFit Receivers



**Hollow Canal** 

RIC

Helix Occluded

Cable Lengths:

1, 2, 3, 4, 5

L, M, P, UP

Receivers:

SnapFit 2.0

SnapFit

Original (two-pin connector)

(ten-pin connector)

(six-pin connector)

Gain:





Canal Lock

Custom Cased or Absolute Power (AP) Receivers Digital SLS

RIC

Custom Cased /

AP Canal







\*Only available in clear

Pro tip: When

measuring the receiver wire for

a Custom Cased

report the number

with the top of the

or AP receiver,

that correlates

entrance to the

ear canal



5,5+

2 3.3+ 4,4+

**BTE Solutions** 

Tubing Lengths:



Thin Tube Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Shell





Skeleton



Canal

Helix

Traditional Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Tubina:

13 Regular

13 Heavy

13 Extra Heavy

13 Dry Tubing 13 Regular Snap

Connect 13 Heavy Snap Connect





Canal

Stent Mold

Digital SLS

Digital SLS

ear canal open.

Telephone Headset



# Earplug Solutions (Hearing Conservation)



#### Solid Earplug

HS Silicone 60 / HS Silicone 40 Affordable solution for

HS Silicone 60 / HS Silicone 40

Contains a patented chamber

to gradually attenuate sound

to appropriate levels (19-22



High-Frequency

Noise Filter\*



RIC

Custom Cased /

AP Canal Lock

#### Swim Plug

Floatable Silicone Custom-fit earplugs that seal out water: floatable if dropped in the water.



#### Sleeping Earplug

HS Silicone 60 HS Silicone 40

Smallest earplugs available using a soft material with a low profile for comfort while sleeping (27 dB NRR).



#### Mini Filtered Earplug

HS Silicone 40

Perfect for musicians, concerts, dental workers, sporting events, noisy work environments, etc. Available in three attenuation levels: 10 dB, 17 dB and 26 dB.

A fourth filter option is our Impulse Filter, providing the wearer with protection from high intensity sounds (IPIL 33.1 dB).

# **Specialty Solutions**



Newscaster/Teleprompter Mold Digital SLS / HS Silicone 60 HS Silicone 40

Skeleton

Used extensively by broadcast professionals, reporters and numerous law enforcement agencies.



#### Racing Receiver Mold

HS Silicone 40

Used by amateur and professional drivers who demand quality one-way listening communication.



#### Pilot Headset with Boom

Digital SLS

Available in two styles: 1. Receiver style earpiece with a snap-on boom; 2. Thick shell style earmold with a screw-set adaptor for the boom.



#### **HA Listening Test Mold** Digital SLS

A hollowed shell of canal

style earmold to hold the

Custom earphone earmolds

and Bluetooth® headsets

are perfect for most cell phone

requiring a mini style adaptor

A shell or skeleton style earmold with a tubing attachment for listening to ITE/ITC/CIC/RIC style and BTE hearing aids for listen checks.



#### Musician's Earplugs HS Silicone 40

**CEP Connection** 

Discreet and cosmetically appealing, these provide attenuation across all frequencies rather than just low and mid-range frequencies.



#### Anesthesiologist Mold Digital SLS

Allows a metal snap-on adaptor to fit 1/8 inch tubing; comes with 42-inch tubing.



#### EarPod/AirPod HS Silicone 60 HS Silicone 40

Earmolds to fit the standard size earbud with a select-avent for acoustic tuning



# **Business Solutions**

In addition to our industry-leading products, we offer a full suite of business solutions designed to keep your practice thriving.



# **Business Equipment**

State-of-the-art diagnostic equipment is essential to providing the best patient experience at your practice. With the help of our Business Equipment Services, it's never been easier to find and maintain the right solutions for all your testing and fitting needs—all in one place. **Equipment includes (but not limited to):** 

- Audiometers
- Tympanometers
- Fitting equipment
- Combination units
- Video otoscopes
- Supplies
- Sound rooms/enclosures

# **All Make Repair**

Our All Make Repair department services hearing aids that are damaged and out of manufacturer's warranty. This even includes hearing aids by manufacturers no longer in business. Our team is here to help with the following:

- All Make Repair estimates
- Reconditioned aids orders
- Miscellaneous parts and accessory orders
- Analog aids
- Body aids
- Bone conduction aids
- All Make Salvage Program
- Hearing Aid Recycle Program

# **Batteries**

Not only does every battery purchase deliver a positive jolt to your bottom line, each visit to your office for batteries is one more chance to build loyalty and generate referrals from your patients.

#### Drive retention

Returning patients are more likely to purchase upgrades, warranties and accessories from you due to the ongoing care and service that you provide.

#### Build loyalty

Superior aftercare will lock in your existing patients rather than allowing them to go to your competition for their future hearing healthcare needs.

#### Boost prospects

Use battery promotions to attract new patients to your practice.

# Protect patients' investment

High quality batteries protect your patient's investment and enhance your patient/provider relationship.

#### Make the right selection

Choose from long-lasting customized or branded batteries, depending on your business goals.



Scan to purchase.

# THE JOY OF HEARING

Audibel logo, Audibel, TeleHear, Pro Fit, Audibel Vitality, Edge Mode, Multiflex Tinnitus, StarLink, Audibel, Neuro Sound Technology, Hear Share, Intrigue, Thrive, SoundGear logo, SoundGear, Phantom, Arc, The Joy of Hearing and Tunz are trademarks of Starkey Laboratories, Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by the Bluetooth SIG. Any use of such marks by Starkey is under license. Other trademarks and trade names are those of their respective owners.

Apple Watch is a trademark of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. "Made for iPod," "Made for iPod," "Made for iPod," "Made for iPod," and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, iPad or Apple Watch may affect wireless performance. Google Play Store is a trademark of Google, LLC.