

Which wireless solution is right for your patient?



PRODUCTS & STYLES	2.4 GHz TECHNOLOGY		900SYNC TECHNOLOGY	NON-WIRELESS
	GENESIS ^{AI}	evolv ^{AI}	Picasso	GENESIS ^{AI}
	<p>micro RIC R RIC RT RIC 312 CIC 312 ITC/HS R ITE R</p> <p><i>with Neuro Sound Technology</i></p>	<p>micro RIC 312 RIC 312 RIC R BTE R BTE 13 PP BTE 13</p> <p><i>with Starkey Sound</i></p>	<p>IIC CIC ITC/HS ITE</p> <p><i>with Acuity Speech Optimization</i></p>	<p>IIC NW</p> <p><i>with Neuro Sound Technology</i></p>
WIRELESS	2.4 GHz ACCESSORIES (StarLink)		SURFLINK	
	<p>Mini Remote Mic Remote Mic Plus Remote Control 2.0 Table Microphone TV Streamer</p>	<p>My Starkey app (for Genesis AI) Thrive Hearing Control App (for Evolv AI) Hear Share App*</p>	<p>Mini Mobile System Remote Microphone 2 Media 2 Remote Control</p>	
FUNCTIONALITY	HEARING PERFORMANCE		HEARING FEATURES	
	<ul style="list-style-type: none"> > 2-Way Audio > Edge Mode > REM Target Match 		<ul style="list-style-type: none"> > Remote Control > TV Streamer > Phone Calls & Audio 	
KEY FEATURES	PERSONAL HEALTH MONITORING		TINNITUS STIMULUS CONTROL	
	<ul style="list-style-type: none"> > Activity Tracking > Engagement Tracking > Fall Detection and Alerts 		<ul style="list-style-type: none"> > Tinnitus Stimulus Control > Noise Control > Remote Microphone 	
KEY FEATURES	INTELLIGENT ASSISTANT		2.4 GHz & 900SYNC TECHNOLOGY	
	<ul style="list-style-type: none"> > Tap Control > Thrive Assistant > Transcribe > Translate > Reminders > Find My Phone 		<ul style="list-style-type: none"> > Industry-best feedback cancellation > Best-in-class adaptive noise management > Ear-to-ear synchronized user controls > Continuously adaptive directionality > Smart frequency-lowering technology > Multiflex Tinnitus Technology > Dedicated music compression focused on listening enjoyment 	
KEY FEATURES			Wireless CROS/BiCROS**	
			<ul style="list-style-type: none"> > REM Target Match 	

*New look. New name. Thrive Care is now Hear Share.

**Available with Genesis AI RIC RT, RIC312 and Evolv AI RIC R, RIC 312, and BTE 13 styles.

2.4 GHz ACCESSORY GUIDE

Bluetooth® Technology
SNR Enhancement and Comfort Features
Personalization Tools
Audio Streaming
Hearing Aid Controls
Convenience Features

Mobile Phone Call Streaming
Media Streaming
Surround Remote Microphone
Focus Remote Microphone
On-Demand Noise Control
User-Created Custom Memories
Automatic Memory Engagement (Geotag)
Noise Control
Automatic Car Memory
Automatic Phone Memory
Multiflex Tinnitus Adjustments
TV Streaming
Line-In Audio Streaming from TV, Computer, MP3, etc.
2-Way Audio
Hearing Aid Volume Changes
Hearing Aid Memory Changes
Mute
Home Shortcut
Multiflex Tinnitus ON/OFF
Find My Hearing Aids
Accessory Volume
Accessory Start/Stop
Teleloop
FM Receiver/DAI

APPLE®	ANDROID™	MINI REMOTE MIC	REMOTE MIC PLUS	REMOTE CONTROL	TABLE MICROPHONE	TV STREAMER
●	● ²		●			
●	● ²		●			
			●		●	
● ³		●	●		●	
●	●			●		
●	●					
●	●					
●	●					
●	●					
		●	●		●	●
			●		●	●
● ⁴						
●	●			●		
●	●			●		
●	●			●		
				●		
●	●			●		
●	●					
●	●					
●	●		●	●		●
●	●			●		
			●			
			●			

● Standard feature ● Feature available with accessory indicated

900SYNC ACCESSORY GUIDE

Bluetooth Technology (iOS and Android Compatible)
SNR Enhancement and Comfort Features
Audio Streaming
Hearing Aid Controls
Convenience Features

Mobile Phone Call Streaming (Apple and Android)
Audio Streaming (Apple and Android)
Surround Microphone
Focus Microphone
On-Demand Noise Control
TV Streaming
Line-In Audio Streaming from TV, Computer, etc.
Volume Changes
Memory Changes
Mute
Home Shortcut
Multiflex Tinnitus ON/OFF

MINI MOBILE SYSTEM	REMOTE MICROPHONE 2	MEDIA 2	REMOTE CONTROL	
			ADVANCED	BASIC
●				
●				
●	●			
●	●			
			●	
●	●	●		
●	●	●		
		●		
			●	●
			●	
			●	
			●	●
			●	

1 Available on Evolv AI standard products except micro RIC 312. | 2 Direct Bluetooth® streaming with Google Pixel 3, 3XL, 3a, 3aXL, 4, 4XL and Samsung Galaxy® S9, S9+, S10 and S10+, Note 10, S20 (requires Android 10). Remote Microphone Plus required for streaming with other Android devices. | 3 Available in Evolv AI to Apple users fit binaurally with low frequency hearing loss ≥ 50 dB. | 4 Available with iPhone 11 and iOS 15 or newer. | Starkey logo, Evolv logo, Thrive logo, TeleHear , Multiflex tinnitus and Edge Mode are registered trademarks of Starkey Laboratories, Inc. Picaso logo, Genesis logo, Neuro Sound Technology, Starkey Sound, StarLink , SurLink logo, Hear Share and Evolv logo are trademarks of Starkey Laboratories, Inc. Apple, the Apple logo, iPhone, iPad, and App Store are trademarks of Apple Inc., registered in the U.S. and other countries. | Android and Google Play are trademarks of Google LLC. | App may have slight differences depending on your phone. | The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Starkey is under license. | iOS is a trademark or registered trademark of Cisco in the U.S. and other countries and is used under license. | ©2023 Starkey Laboratories, Inc. All Rights Reserved. 81433-600 8/23 FLYR2728-18-EN-GS

