

Auracast assistant



Connect anywhere with **Auracast assistant**

Stay connected wherever you go! Easily find and join Auracast™ broadcasts—whether you're at the airport, the gym, or listening to public announcements. Enjoy instant, clear audio streamed directly to your hearing aids from the My Audibel app.

ARIS^{AI}
VITALITY^{AI}



Where to find it

My Audibel mobile app*

*Requires Hearing Aid Firmware version 10.3.1.5 or later.

Connecting to Auracast broadcasts through the My Audibel mobile app

1. From the Home screen, tap the **accessories icon**.
2. Tap **Auracast™ broadcasts**.
3. Tap the broadcast you would like to join.
4. Experience seamless audio streaming to your compatible Audibel hearing aids.

*Scan to check
device compatibility*

