

Wireless **Accessories**

Discover our innovative lineup




A connected experience is an **expectation**.

Wireless accessories stream the sounds you love directly to your hearing aids.

Our accessories, paired with our wireless hearing aids, connect you to your favorite music, TV shows, and phone calls with superior sound quality.

Enabled with Bluetooth® technology, they are easy to pair and even easier to use.





StarLink Edge Edge TV Streamer

Future-proof and *Auracast*™ ready, the Edge TV Streamer uses
Bluetooth® Low Energy (LE) Audio protocol to stream audio straight from a TV to paired LE Audio-enabled hearing aids, like Vitality AI.




StarLink Remote Control 2.0

Redesigned for maximum simplicity, the remote lets users easily control volume, change programs, and quickly access a favorite feature. Includes replaceable cell battery and lanyard.

Easily connect to your life.

Talk to your hearing care professional about finding the right accessories for your most important hearing needs.





StarLink Remote Microphone+

In noisy environments, easily stream to both ears from Bluetooth, Telecoil, FM or Line-In and enjoy your one-on-one conversations to the fullest.



StarLink TV Streamer

Easily stream audio from your TV or other electronic audio source directly to your hearing aids. It offers excellent sound quality, is simple to use, and supports both analog and digital input sources.



StarLink Table Microphone

A multi-functional device designed to improve listening in group gatherings or noisy restaurants. It also works as a TV streamer or even a body-worn remote microphone.



StarLink Mini Remote Microphone

Enjoy one-on-one conversations in noisy environments with our small, easy-to-use Mini Remote Microphone. Simply clip it to the clothing of the person you're talking to and enjoy the conversation.

StarLink, Audibel and the Audibel logo are registered trademarks of Starkey Laboratories, Inc. Audibel Vitality is a trademark of Starkey Laboratories, Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by the Bluetooth SIG. and any use of such marks is under license. Other trademarks and trade names are those of their respective owners.



www.Audibel.com

6700 Washington Ave. S. Eden Prairie, MN 55344