

ALL NEW

ARISAI

Technology so intelligent, it's superhuman



"I am so confident and completely aware of everything."

Bill P. *Photographer*

Every patient deserves to hear their best, in a way that feels natural and true. You strive to provide hearing aid wearers like Bill with a solution that allows them to let go of their hearing loss insecurities and take back their lives—navigating every situation with a genuine sense of connection and clarity.

All-new Aris Al is that solution.



Watch Bill's story



Introducing DNN360

First, Audibel introduced On-Demand DNN. Then, Always-On DNN. Now, with Aris AI, we present the world's first DNN-powered directionality system, DNN 360.

Aris Al intuitively detects characteristics such as speech, background noise, and motion, adjusting accordingly. This helps patients navigate their surroundings with greater confidence, knowing their hearing will adapt as they go.



Omnidirectional



ALL NEW

DNN Spatial

Awareness

Aris Al automatically applies the most appropriate directionality patterns for the best listening experience in the most challenging situations—without compromising on spatial awareness.



ALL NEW
DNN Directionality



Motion Based Optimization



"I am so confident and completely aware of everything."


Bill P. *Photographer*

Every patient deserves to hear their best, in a way that feels natural and true. You strive to provide hearing aid wearers like Bill with a solution that allows them to let go of their hearing loss insecurities and take back their lives—navigating every situation with a genuine sense of connection and clarity.

All-new Aris Al is that solution.



Watch Bill's story



Leading the way in AI, we offer Gen AI features in the My Audibel app and Pro Fit software that open new paths to support you and your patients.

Gen Al Smart Assistant


Gen AI Smart Assistant allows patients to get **hands-free help** using their voice. They can receive detailed answers to questions, switch programs, set reminders, and more.

TeleHear AI INDUSTRY FIRST

A world first—TeleHear AI allows hearing aid wearers to manage their own sound concerns. After listening to the patient's issue, the app feature **instantly creates a sound adjustment** for them, based on professional expertise.

Al Assist ALL NEW

Al Assist answers your product-related questions by leveraging Gen Al and current Audibel content.



We are continuously innovating new features that can help you improve every patient's hearing journey and quality of life.

Hearing health is connected to overall well-being, which is why we continue to pioneer new health and wellness features available to patients using the My Audibel app.

Wellness at their fingertips

Balance Builder ALL NEW

Balance Builder is an **all-new feature** that helps patients proactively prevent falls through physical awareness exercises designed to bolster balance and stability.

Respiratory Rate monitoring INDUSTRY FIRST

Respiratory Rate, an **industry-first feature**, allows patients to monitor their own respiratory rate—one of the four vital signs that help measure the body's most basic functions.



"I am so confident and completely aware of everything."


Bill P. *Photographei*

Every patient deserves to hear their best, in a way that feels natural and true. You strive to provide hearing aid wearers like Bill with a solution that allows them to let go of their hearing loss insecurities and take back their lives—navigating every situation with a genuine sense of connection and clarity.

All-new Aris Al is that solution.



Watch Bill's story



Leading the way in AI, we offer Gen AI features in the My Audibel app and Pro Fit software that open new paths to support you and your patients.

Gen Al Smart Assistant


Gen Al Smart Assistant allows patients to get **hands-free help** using their voice. They can receive detailed answers to questions, switch programs, set reminders, and more.

TeleHear AI INDUSTRY FIRST

A world first—TeleHear AI allows hearing aid wearers to manage their own sound concerns. After listening to the patient's issue, the app feature **instantly creates a sound adjustment** for them, based on professional expertise.

Al Assist ALL NEW

Al Assist answers your product-related questions by leveraging Gen Al and current Audibel content.



We are continuously innovating new features that can help you improve every patient's hearing journey and quality of life.

Hear with clarity, connect with confidence



Taking better hearing in a whole new direction, we've engineered Aris Al with our **industry-leading Al and G3 Gen Al Neuro Processor.** The result: innovation designed to come closer to restoring the normal sense of hearing than any hearing technology before it.





Faster fittings, greater satisfaction

Now, with **Acoustic Model Optimization 2.0**, Pro Fit is faster than ever before. AMO 2.0 leverages the feedback cancellation initialization signal and in-situ measurements to estimate the acoustic-coupling characteristics when the hearing aid or receiver is placed in the ear.

When used with Best Fit, AMO 2.0 improves first-fit accuracy and reduces the manual adjustments needed to match targets.



Learn more about Aris Al at AudibelPro.com





Audibel, TeleHear, Edge Mode, Aris, DNN 360, Balance Builder, and the Audibel logo are trademarks of Starkey Laboratories, Inc. The Auracast™ word mark and logos are trademarks owned by Bluetooth SIG, Inc. and any use of such marks is under license. ⊚2025 Starkey Laboratories, Inc. All Rights Reserved. P00012821 8/25 BKLT3268-00-EN-AB

^{*}Marquardt, D., Xiao, J., Ganeshkumar, K., Xu, J. Taylor, L., and McKinney, M. (2025) Enhancing Directionality with Deep Neural Networks