

evolv^{AI}



Next level sound. Next level everything.

Evolv AI brings a full line of styles to the world's most innovative hearing aids.

Experience Starkey's award-winning, most successful hearing aid line to date. Evolv AI seamlessly connects and adapts with your patients' lives to deliver effortless hearing in even more listening environments. Now featuring enhanced Edge Mode, improved 2-Way Audio, and health monitoring.

Starkey Sound is effortless.

Performing 55 million automatic adjustments per hour, Evolv AI is reimagining effortless hearing.

Our always-on, always automatic, and always effortless hearing technology delivers realistic and genuine sound quality in every listening environment without the need to do anything extra. Evolv AI now expands its 40% additional reduction in noise energy to include wind and machine noise. Enhanced Edge Mode offers even clearer sound in cars and other vehicles.

	Premium 2400	Advanced 2000	Select 1600	Basic 1200	Entry 1000
Thrive 2.4 GHz + NFMI* Technology					
Sound Imaging: Channels Bands	24	20	16	12	10
Edge Mode	●	◐	◐		
Voice AI	●				
Speech Optimization	●	●	●	●	●
Music Optimization					
Music Adaptation	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
E2E Music Adaptation	●	●	●		
Ear-to-Ear Technology					
E2E Wind Noise Management	●	●	●		
E2E Machine Noise Adaptation	●	●	●		
E2E Directionality**	●	●	●		
E2E Phone Streaming*	●	●	●		
Sound Manager					
Auto Music	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Quiet	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4	1 2 3 4
Spatial Speech Enhancement	●	●	●		
Machine Noise	Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction	Up to 7 dB of reduction
Speech in Noise	Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction
Transient Noise Reduction	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 8 dB of reduction	Up to 8 dB of reduction	Up to 6 dB of reduction
Wind	Up to 32 dB of reduction	Up to 17 dB of reduction	Up to 9 dB of reduction	Up to 9 dB of reduction	Up to 9 dB of reduction
Directional Processing					
Immersion	●	●	●	●	
Adaptive	●	●	●	●	
Dynamic	●	●	●	●	●
Telehealth Service					
TeleHear	●	●	●	●	●
Telecoil					
Telecoil	Available on ITE R, RIC R, BTE R and BTE 13.				

*Available on RIC R, RIC 312 and BTE 13. **Available for products with directional microphones.

A complete style portfolio makes effortless hearing the expectation.


We're combining our exceptional Starkey Sound with the world's most innovative complete family of hearing technology.

Custom craftsmanship and innovation have led to an effortless hearing experience. Patients now have a complete line of hearing aid options that fit seamlessly into their unique lifestyles. From our smallest to our most powerful, we've cracked the code for our hearing aids to be effortlessly perfect.



A range of styles for your patients' specific needs.

Effortlessly fit every style for every patient.



IIC NW, CIC NW, CIC, ITC/HS R, ITE R, micro RIC 312, RIC 312, RIC R, BTE R, BTE 13, Power Plus BTE 13

Clearer sound. Less background noise. Effortless hearing.

Evolv AI introduces our most advanced hearing aids that seamlessly fit into your patients' lives.

- 2-Way Audio** ENHANCED
Included in every wireless style, 2-Way Audio now allows patients to take and end calls directly from their Evolv AI hearing aids instead of their Apple device.
- Battery Performance** ENHANCED
Re-engineered for 17% higher rechargeable battery capacity and up to 50% less current drain in noisy conditions.

Edge Mode ENHANCED
Putting AI into every patient's hands with on-demand features like Edge Mode has proven speech audibility advantages. Now enhanced to offer even clearer sound in cars and other vehicles.

Fall Alert
With significant improvements and expanded datasets, Fall Alert now provides even more accuracy in detection than ever before.

Reminders
Patients get reminder and calendar alerts directly through their hearing aids.

TeleHear ENHANCED
Easier-to-use remote programming offers even more flexibility and fine tuning when virtually making adjustments to patients' hearing aids.

Do more with Thrive
With the industry's first health monitoring feature
The Thrive Hearing Control app delivers powerful features that are designed to help patients monitor activities to live healthier — while Hear Share provides peace of mind for caregivers.



*Compatible with iPhone 11, iPhone 11 Pro, iPhone 11 Pro Max, and later, as well as iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), and iPad mini (6th generation).

Effortless intelligence

Our latest technology provides superior sound quality plus advanced features and capabilities.

	Premium 2400	Advanced 2000	Select 1600	Basic 1200	Entry 1000
Voice AI	●				
Translate	●				
Multiflex Tinnitus Pro	●				
Edge Mode: On-demand control	●				
Transcribe	●	●			
Edge Mode: Thrive app access	●	●	●		
Thrive Assistant	●	●	●		
Find My Phone	●	●	●		
Engagement & Activity Tracking	●	●	●	●	
Fall Alert	●	●	●	●	
Hear Share app compatibility	●	●	●	●	
CROS system*	●	●	●	●	
2-Way Audio	●	●	●	●	●
Reminders	●	●	●	●	●
Self Check	●	●	●	●	●
Auto On/Off	●	●	●	●	●
Tap Control	●	●	●	●	●
TeleHear full remote fit	●	●	●	●	●
REM Target Match	●	●	●	●	●

Accessories your patients love

Patients connect effortlessly to their world with our versatile line of accessories.



Table Microphone
Make group settings easy with eight built-in microphones that stream the primary speaker's voice directly to your patient's hearing aids.



Remote Microphone + & Mini Remote Microphone
Enables easy audio streaming and one-on-one conversations in noisy environments.



TV Streamer
Allows patients to easily stream audio from their TV or other electronic audio source directly to their hearing aids.



Remote
Easily control memory and volume, mute hearing aids and turn other special features on and off.



Standard & Custom Charger
With up to three full charges, our chargers let patients know when their rechargeable hearing aids are charged, automatically turning the devices on when removed from the charger.



Mini Turbo Charger
In just seven minutes, our pocket-sized charger provides 3.5 hours of power and can completely charge RIC R or BTE R hearing aids up to four times.

*Available on RIC R, RIC 312 and BTE 13.

**Learn how you can continue to help
your patients live their lives to the
fullest at starkey.com.**



StarkeyPro.com