



Product Portfolio Guide

Spring 2026



Tech that's yesterday's thinking **or future driven?**

From enhancements to our in-demand Omega AI to a new essential line featuring two new BTEs, Starkey is setting the speed of innovation with technology that lets you fit AI-powered performance on every patient—all using the same fitting software and app.

ASK MORE from your technology partner

ENHANCED
OMEGA^{AI}

Tech Tier/Channels/Bands

24 / 20 / 16

AURACAST

CROS COMPATIBLE



mRIC R
Waterproof, Rechargeable, LED, Wireless

RIC RT
Waterproof, Rechargeable, LED, Wireless, Telecoil

RIC 312
Wireless



CIC
Wireless

ITC R
Waterproof, Rechargeable, Wireless

ITE R
Waterproof, Rechargeable, Wireless, Telecoil



- Waterproof*
- Rechargeable
- LED Status Indicator Light
- Wireless
- Telecoil

*Up to 1m.

Available Colors

- Standard Colors**
- White
 - Tech Black
 - Graphite Gray
 - Silver
 - Chestnut
 - Caramel
 - Beige


- Custom Colors**
- Faceplate Colors**
- Pink
 - Light Brown
 - Medium Brown
 - Chestnut
 - Dark Brown
 - Black

- Shell Colors**
- Clear
 - Pink
 - Light Brown
 - Blue/Red
 - Black

SnapFit 2.0 Stock and Custom Cased Receivers

Comfort Bud Open 5 mm, 6 mm, 7 mm, 9 mm, 11 mm	Comfort Bud Occluded 4 mm, 5 mm, 6 mm, 7 mm, 9 mm, 11 mm	Power Dome 8-10 mm, 10-12 mm, 12-14 mm	Custom Cased Receiver Custom shell with embedded receiver
Stock Silicone Helix RIC Earmold	Stock Hollow Silicone RIC Earmold	Stock Silicone RIC Earmold	

Fitting Ranges



Visit StarkeyPro.com/OmegaAI for full technical data sheets.



Engineered with our industry-leading AI, our G3 Gen AI Neuro Processor, and the world's first DNN-powered directionality system, DNN 360, Omega AI is Starkey's most advanced hearing technology yet. It's designed to come closer to restoring the normal sense of hearing than any hearing technology before it.

For more Omega AI details, visit StarkeyPro.com/OmegaAI



OMEGA^{AI}

Tech Tiers | Channels | Bands

24

20

16

Edge Mode+				
Automatic & User Control Personalization <i>App Control</i>		● ● ●	○ ● ●	○ ○ ●
Speech Optimization				
Consonant Brightness <i>Voice AI*</i>		● ● ● ● ● ●	● ○ ○ ○ ○ ○	● ○ ○ ○ ○ ○
Sound Manager				
DNN 360	DNN Directionality**	●	○	○
	DNN Spatial Awareness**	●	○	○
	Quiet	Maximum expansion ●●●●	Moderate expansion ●●○	Mild expansion ●○○
	Machine Noise Reduction	Up to 22 dB of reduction ●●●●	Up to 12 dB of reduction ●●○○	Up to 7 dB of reduction ●○○○
	Speech in Noise Reduction	Up to 22 dB of reduction ●●●●	Up to 10 dB of reduction ●●○○	Up to 8 dB of reduction ●○○○
	Transient Noise Reduction	Up to 15 dB of reduction ●●●●	Up to 9 dB of reduction ●●○○	Up to 6 dB of reduction ●○○○
	Wind Reduction	Up to 35 dB of reduction ●●●●	Up to 19 dB of reduction ●●○○	Up to 10 dB of reduction ●○○○
Additional Key Features				
Gen AI Smart Assistant & TeleHear AI Translate		●	○	○
<i>ALL NEW</i> Push-to-talk		●	○	○
Transcribe		●	●	○
Multiflex Tinnitus Pro		●	●	○
<i>ALL NEW</i> Auracast™ assistant		●	●	●
<i>ALL NEW</i> Google Fast Pair™		●	●	●
Smart Assistant & 2-Way Audio		●	●	●
Hearing Health & Activity Tracking		●	●	●
Fall Alert & Balance Assessment		●	●	●
Reminders & Self Check		●	●	●
Tap Control		●	●	●
TeleHear full remote fit		●	●	●
REM Target Match		●	●	●
Demo Stream		●	●	●

Wireless Accessories



StarLink Edge LE Adapter
Part No. P00018678



StarLink Edge Remote Microphone
Part No. P00004321



StarLink Edge Table Microphone
Part No. P00004326



StarLink Edge TV Streamer
Part No. P00004315



StarLink Remote Control 2.0
Part No. P00002347



Charging Solutions

All chargers include Turbo Charge

Serve as a charging station and carrying case for your rechargeable hearing aids.



StarLink Premium Charger 2.0 w/ onboard battery

- mRIC R
- RIC RT



StarLink Premium Custom Charger 2.0 w/ onboard battery

- ITC/ITE R


StarLink Charger 2.0

- mRIC R
- RIC RT


StarLink Custom Charger 2.0

- ITC/ITE R

Compatible Apps



My Starkey



Hear Share



Scan QR code to check phone compatibility



Tech Tier/Channels/Bands
24 / 20* / 16*



Waterproof**
Rechargeable

For more product details,
visit StarkeyPro.com



Faceplate Colors

- Pink
- Light Brown
- Medium Brown
- Chestnut*
- Dark Brown
- Black

Shell Colors

(IIC only available in Black)

- Clear
- Pink
- Light Brown
- Blue/Red
- Black

Charging Solutions


Both chargers include
Turbo Charge



StarLink Premium Custom Charger 2.0
w/ onboard battery
• CIC R NW

StarLink Custom Charger 2.0
• CIC R NW

Fitting Ranges



Tech Tiers Channels Bands	24	20*	16*
Edge Mode+ (CIC styles only)			
User Control	●	○	○
Speech Optimization			
Consonant Brightness	● Maximum	● Mild	● Mild
Sound Manager			
Quiet	Maximum expansion	Moderate expansion	Mild expansion
Machine Noise Reduction	Up to 22 dB of reduction	Up to 12 dB of reduction	Up to 7 dB of reduction
Speech in Noise Reduction	Up to 22 dB of reduction	Up to 10 dB of reduction	Up to 8 dB of reduction
Transient Noise Reduction	Up to 15 dB of reduction	Up to 9 dB of reduction	Up to 6 dB of reduction
Wind Reduction	Up to 35 dB of reduction	Up to 19 dB of reduction	Up to 10 dB of reduction
Additional Key Features			
Multiflex Tinnitus Pro with new tinnitus maskers	●	●	○
REM Target Match	●	●	●
Waterproof**	●	●	●



Powered by our industry-leading Neuro Sound Technology, Starkey Signature Series delivers the size, performance and personalized fit that patients can only get from you.

*CIC styles only. **Up to 1m. Rechargeable style only.

Quality in... quality out

Your professional expertise is critical when selecting the best solution for a patient. Many factors and considerations come into play when determining if a patient is a custom hearing aid candidate, including:

- Audiogram
- Dexterity
- Ear anatomy and history
- Lifestyle

Ultimately, **the patient has come to you** to weigh all considerations—physical, audiological and lifestyle—and recommend a solution that delivers both success and satisfaction.

A superior custom product fitting starts with a high quality ear impression. Making a perfect impression requires proper technique, attention to detail, and repetition.

As the custom hearing aid experts, we know a thing or two about making accurate, precision-tested impressions—and how good impressions can also be used to make earmolds, custom cases, hearing protection solutions and more.



Scan the QR code or go to StarkeyPro.com for support materials, how-to videos, and more.

Developed with help from Starkey's most experienced impression takers, these tools provide best practice tips perfected from shooting thousands of impressions in our Center for Excellence.

EDGE^{AI}

Tech Tier/Channels/Bands

24 / 20 / 16

AURACAST

CROS COMPATIBLE



mRIC R



RIC RT



RIC 312



- Waterproof*
- Rechargeable
- Wireless
- Telecoil

*Up to 1m.

Available Colors

- Standard Colors
- White
 - Tech Black
 - Graphite Gray
 - Silver
 - Chestnut
 - Caramel
 - Beige

SnapFit 2.0 Stock and Custom Cased Receivers



Comfort Bud Open
5 mm, 6 mm, 7 mm,
9 mm, 11 mm



Comfort Bud Occluded
4 mm, 5 mm, 6 mm,
7 mm, 9 mm, 11 mm



Power Dome
8-10 mm, 10-12 mm,
12-14 mm



Custom Cased Receiver
Custom shell with
embedded receiver



**Stock Silicone Helix
RIC Earmold**





**Stock Hollow
Silicone RIC Earmold**



**Stock Silicone
RIC Earmold**

Fitting Ranges



Visit StarkeyPro.com/EdgeAI for full technical data sheets.



Powered by the G2 Neuro Processor with the neural processing unit (NPU) integrated into the chip, and Bluetooth® Low Energy (LE) Audio streaming protocol, Starkey's Edge AI advanced processing is able to provide so much more than all-day battery life and waterproof protection – allowing patients to live easier, more connected lives.

For more Edge AI details,
visit StarkeyPro.com/EdgeAI























EDGE^{AI}

Tech Tiers | Channels | Bands

24

20

16

Edge Mode+				
Automatic	●	○	○	○
User Control	●	○	○	○
Personalization	●	●	○	○
App Control	●	●	●	●
Speech Optimization				
Consonant Brightness	●	●	●	●
Voice AI*	●	○	○	○
Sound Manager				
Quiet 	 Maximum expansion	 Moderate expansion	 Mild expansion	
Machine Noise Reduction 	 Up to 22 dB of reduction	 Up to 12 dB of reduction	 Up to 7 dB of reduction	
Speech in Noise Reduction 	 Up to 22 dB of reduction	 Up to 10 dB of reduction	 Up to 8 dB of reduction	
Transient Noise Reduction 	 Up to 15 dB of reduction	 Up to 9 dB of reduction	 Up to 6 dB of reduction	
Wind Reduction 	 Up to 35 dB of reduction	 Up to 19 dB of reduction	 Up to 10 dB of reduction	
Additional Key Features				
Translate	●	○	○	○
Transcribe	●	○	○	○
Multiflex Tinnitus Pro	●	●	○	○
ALL NEW Auracast™ assistant	●	●	●	●
Smart Assistant & Find My Hearing Aids	●	●	●	●
2-Way Audio	●	●	●	●
Hearing Health & Activity Tracking	●	●	●	●
Fall Alert & Balance Assessment	●	●	●	●
Reminders & Self Check	●	●	●	●
Tap Control	●	●	●	●
TeleHear full remote fit	●	●	●	●
REM Target Match	●	●	●	●
Demo Stream	●	●	●	●
Waterproof up to 1m**	●	●	●	●

Wireless Accessories



StarLink Edge LE Adapter
Part No. P00018678



StarLink Edge Remote Microphone
Part No. P00004321



StarLink Edge Table Microphone
Part No. P00004326



StarLink Edge TV Streamer
Part No. P00004315



StarLink Remote Control 2.0
Part No. P00002347



Charging Solutions

All chargers include Turbo Charge



StarLink Premium Charger 2.0 w/ onboard battery


- mRIC R
- RIC RT




StarLink Charger 2.0

- mRIC R
- RIC RT

Compatible Apps



My Starkey



Hear Share



Scan QR code to check phone compatibility



INTRODUCING

G SERIES AI

Tech Tier/Channels/Bands

24* / 20* / 16* / 12 / 10



- Waterproof**
- Rechargeable
- Wireless
- Telecoil


*BTE styles only.
**Up to 1m.

Available Colors





Standard Colors

- Graphite Gray
- Silver
- Chestnut
- Beige

SnapFit 2.0 Stock and Custom Cased Receivers



Fitting Ranges



Visit StarkeyPro.com for full technical data sheets.



Powered by our industry-leading Neuro Sound Technology, G Series AI offers an easy, affordable, and immediate way to improve patient hearing and reconnect them with the world around them.







For more product details, visit StarkeyPro.com



G SERIES^{AI}

Tech Tiers | Channels | Bands

24* **20*** **16*** **12** **10**

Speech Optimization					
Ear-to-Ear Technology	●	●	○	○	○
Frequency Lowering	●	●	●	●	●
Sound Manager					
Directionality 	Adaptive	Adaptive	Adaptive	Dynamic	Dynamic
Transient Noise 	●	○	○	○	○
Quiet 	● Maximum expansion	●	●	●	●
Machine Noise Reduction 	● Up to 22 dB of reduction	● Up to 7 dB of reduction	● Up to 7 dB of reduction	● Up to 7 dB of reduction	● Up to 7 dB of reduction
Speech in Noise Reduction 	● Up to 22 dB of reduction	● Up to 10 dB of reduction	● Up to 8 dB of reduction	● Up to 8 dB of reduction	● Up to 8 dB of reduction
Wind Reduction 	● Up to 35 dB of reduction	● Up to 19 dB of reduction	● Up to 10 dB of reduction	● Up to 10 dB of reduction	● Up to 10 dB of reduction
Additional Key Features					
Edge Mode+	●	●	○	○	○
Reminders	●	●	○	○	○
Multiflex Tinnitus Technology	●	● (white, pink, & ocean noise only)	● (white, pink, & ocean noise only)	● (white noise only)	○
Find My Hearing Aids	●	●	●	●	●
Self Check	●	●	●	●	●
TeleHear full remote fit	●	●	●	●	●
REM Target Match	●	●	●	●	●
Apple Watch Controls**	●	●	●	●	●

Wireless Accessories




Charging Solutions

All chargers include
Turbo Charge



Compatible Apps



My Starkey



Scan QR code to check
phone compatibility



Features noted above in blue are available within My Starkey app.

*BTE styles only. **For iOS users only.

My Starkey App

Easy as it gets. Powerful as you need.



MyStarkey

Focus on benefits

Omega AI, Edge AI, and G Series AI wireless hearing aids can connect via Bluetooth® technology to the My Starkey app.

- **Hear better anywhere**

Make instant adjustments for clearer speech and reduced noise – even in the toughest listening environments.

- **Get tips, tools, and resources**

Adapt to hearing aids faster with quick how-to's, tips, and videos. All tailored for where the hearing aid wearer is on their hearing journey.




- **Monitor health and wellness***

Monitor daily hearing and physical activity to help stay healthy. Set and achieve goals for hearing, steps, exercise, and beyond.


- **All from the palm of your hand or wrist**

Change programs and adjust volume quickly and easily on a compatible Apple Watch. Visit the Learn section within My Starkey for installation instructions.

How to download

1. Open the **App Store**  on Apple devices or **Google Play**  on Android devices.
2. In the search box, type  **My Starkey**.
3. Tap **GET** on Apple devices or **INSTALL** on Android devices. This will download the app to the compatible device.





A new dimension in connectivity

All-new **Auracast™ assistant** and **Google Fast Pair™** technology make connectivity so convenient, it's automatic.

ALL NEW

To learn more, visit StarkeyPro.com



*Compatible with Omega AI and Edge AI products only.

Other Apps


Enhance the capabilities of our hearing aids and simplify daily life with our lineup of popular hearing apps.



SoundCheck Live

Helps anyone gain awareness of noise exposure and the importance of protecting hearing by:

- Measuring sound levels around them using the built-in sound level meter.
- Learning about environmental noise's impact on hearing and how hearing protection devices like SoundGear can help.



Hear Share

Works with My Starkey app to allow wearers of Omega AI and Edge AI to:

- Share information about hearing aid usage, social and physical activity progress, safety, and more—with the people they choose.
- Provide peace of mind for patients and their caregivers.



Hearing Protection

Be your community's full-service hearing care expert by offering hearing protection products from SoundGear. SoundGear products come in a variety of electronic and non-electronic solutions.



Electronic Solutions



PHANTOM Custom Fit Electronic

World's first custom, rechargeable, and Bluetooth® compatible hearing aid and protection device. Up to 30 dB gain.



SHIELD

Small, light and ready to take everywhere delivering dynamic digital hearing protection anytime you need it.






 **93 dB** | Part No. P00001867

 **80 dB** | Part No. P00001866



RACING RECEIVERS

Custom design allows for maximum comfort and durability for amateur and professional drivers to experience quality one-way listening communication.

	Top Seller	
	Phantom	Shield
Fit	Custom	Universal
Rechargeable	Yes	No
NRR	22 dB	26 dB
Warranty	1 year	1 year
Color Options		
Lanyard		
Bluetooth		
Digital Modes		

Non-electronic Solutions



Filtered Earplug HS Silicone

Referred to as musician or industrial earplugs, these earplugs attenuate across all frequencies.



Solid Earplug HS Silicone

Outstanding noise protection in a variety of recreational and occupational environments.



High-Frequency Noise Filter HS Silicone

Referred to as "Hocks", these earplugs contain a patented acoustical chamber and filter that attenuate sound to appropriate levels.

Additional SoundGear earmolds:

Newscaster/Teleprompter mold
Digital SLS/HS Silicone

Bluetooth mold
HS Silicone

Pilot headset with boom
Digital SLS

Sleeping earplug
HS Silicone

Anesthesiologist mold
Digital SLS

Swim earplug
Floatable Silicone

For a full selection of custom options, please contact Customer Service.

Personal Audio

For musicians, music aficionados, or professionals looking for a custom audio solution, Starkey has options to maximize their hearing experience.

Tunz Duo S | Tunz Duo | Tunz Solo

Additional Items

SoundGear Custom Charger for Phantom devices
Part No. 43059-610

Replacement Tri-flanges for Shield
Size 1: P00009811 (4pk)
Size 2: P00009812 (4pk)

Foam Buds Assorted for Shield

Part No. P00009810 (6pk)

Replacement Filters

Filtered earplug: 10 dB, 17 dB, 26 dB, Impulse, Total block
High-Frequency (Hock): Part No. 50293-000

SoundGear Black Lanyard

Variety of lengths available

SoundGear Zipper Pouch

Part No. 91289-007

Hearing Aid Earmolds

As the leader in customization, Starkey understands offering your patients an expansive portfolio of hearing care solutions specific to their needs is key.



Hear better. Live better.

Material Options

Digital SLS Acrylic	• RIC • BTE & Thin Tube • Specialty*	• One color	Most Rigid
			90 Shore
High Strength Silicone 60	• RIC • BTE & Thin Tube • Specialty*	• Swirl - max. 3 colors • Glitter option	60 Shore
High Strength Silicone 40	• RIC • BTE & Thin Tube • Specialty*	• Swirl - max. 3 colors • Glitter option	40 Shore
Floatable Silicone	• Swim Earplug (SoundGear only)	• Swirl - max. 3 colors	25 Shore
			Softest

Venting Guide

This guide is applicable for all hearing aid earmolds.** Vent size is limited by ear anatomy.

Vent Name	Approximate Size
Small	1.2 mm
Medium	1.8 mm
Large	2.55 mm
Select-A-Vent	3.37 mm
Extra Large	As large as canal allows

Specialty Solutions



HA Listening Test Mold
Digital SLS

A shell or skeleton style earmold with a tubing attachment for listening to ITE/ITC/CIC/RIC style and BTE hearing aids for listen checks.




Stent Mold
Digital SLS

A hollowed shell of canal style earmold to hold the ear canal open.

RIC Solutions

RIC Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Available for both SnapFit 2.0 and SnapFit Receivers




Custom Cased or Absolute Power (AP) Receivers Digital SLS

Cable Lengths:
1, 2, 3, 4, 5

Gain:
L, M, P, UP

Receivers:
SnapFit 2.0 (ten-pin connector)
SnapFit (six-pin connector)
Original (two-pin connector)




*Only available in clear
Pro tip: When measuring the receiver wire for a Custom Cased or AP receiver, report the number that correlates with the top of the entrance to the ear canal.

Color Options*

Our earmolds come in a variety of color options. Some colors appear differently in certain color families based on material type. **Swirl combinations are also available using up to three colors.**

Digital SLS



HS Silicone 60



HS Silicone 40



Glitter




Pro Tip: Adding White as one of the three colors within a color swirl combination will help separate and make the other color choices "pop."
Color combinations to avoid:
• 'Light' colors combined in Swirl or 'Light' colors combined with 'Clear'
• Similar colors combined in Swirl

BTE Solutions

Thin Tube Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Tubing Lengths:


- 1
- 2
- 3, 3+
- 4, 4+
- 5, 5+



Traditional Earmolds Digital SLS / HS Silicone 60 / HS Silicone 40

Tubing:

- 13 Regular
- 13 Heavy
- 13 Extra Heavy
- 13 Dry Tubing
- 13 Regular Snap Connect
- 13 Heavy Snap Connect



Business Solutions

In addition to our industry-leading products, we offer a full suite of business solutions designed to keep your practice thriving.



Hear better. Live better.

Business Equipment

State-of-the-art diagnostic equipment is essential to providing the best patient experience at your practice. With the help of our Business Equipment Services, it's never been easier to find and maintain the right solutions for all your testing and fitting needs— all in one place. **Equipment includes (but not limited to):**

- **Audiometers**
- **Tympanometers**
- **Fitting equipment**
- **Combination units**
- **Video otoscopes**
- **Supplies**
- **Sound rooms/enclosures**

All Make Repair

Our All Make Repair department services hearing aids that are damaged and out of manufacturer's warranty. This even includes hearing aids by manufacturers no longer in business. **Our team is here to help with the following:**

- **All Make Repair estimates**
- **Reconditioned aids orders**
- **Miscellaneous parts and accessory orders**
- **Analog aids**
- **Body aids**
- **Bone conduction aids**
- **All Make Salvage Program**
- **Hearing Aid Recycle Program**

Batteries

Not only does every battery purchase deliver a positive jolt to your bottom line, each visit to your office for batteries is one more chance to build loyalty and generate referrals from your patients.

- **Drive retention**
Returning patients are more likely to purchase upgrades, warranties and accessories from you due to the ongoing care and service that you provide.
- **Build loyalty**
Superior aftercare will lock in your existing patients rather than allowing them to go to your competition for their future hearing healthcare needs.
- **Boost prospects**
Use battery promotions to attract new patients to your practice.
- **Protect patients' investment**
High quality batteries protect your patient's investment and enhance your patient/provider relationship.
- **Make the right selection**
Choose from long-lasting customized or branded batteries, depending on your business goals.

Hearing is Our Concern

Manufacturer according to Health Canada: Starkey Labs Canada Co. 2476 Argenta Road Suite 301, Mississauga, ON L5N 6M1; www.starkeycanada.ca.

Starkey logo, Starkey, TeleHear, Starkey Edge, Starkey G Series, Edge Mode, Multiflex Tinnitus, StarLink, Starkey Signature Series, SoundGear Shield, Omega, Sound [+checkmark], Hear Share, SoundGear logo, SoundGear, Phantom, Hearing is Our Concern, Hear better. Live better., Neuro Processor, DNN 360 and Tunz are trademarks of Starkey Laboratories, Inc.

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The Auracast™ word mark and logos are trademarks owned by the Bluetooth SIG. Any use of such marks by Starkey is under license. Other trademarks and trade names are those of their respective owners.

Apple Watch and App Store are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. "Made for iPod," "Made for iPhone" and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPod, iPhone, or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPod, iPhone, iPad or Apple Watch may affect wireless performance. Google Play Store, Android and Fast Pair are trademarks of Google, LLC.