

Hear from all directions  
and stay focused while  
actively participating in  
conversations in every  
listening environment.



[www.Audibel.com](http://www.Audibel.com)

Global Headquarters  
6700 Washington Ave. S.,  
Eden Prairie, MN 55344  
1.800.769.2590

Audibel, Audibel logo, Arc and Arc logo are trademarks of Starkey Laboratories, Inc.

©2021 Starkey Laboratories, Inc. All Rights Reserved. 86582-010 8/21 BROCC2942-00-EE-AB

# Arc<sup>AI</sup>

## CROS/BiCROS System

Made for single-sided  
hearing loss



# Single-sided hearing loss solutions are effortless with Arc AI.

**Our CROS/BiCROS System paired with our most advanced hearing technology and Audibel Sound™ is making hearing effortless.**

Arc AI product family offers a wireless CROS/BiCROS system in select styles. Loaded with our most innovative technology and clearest sound, these hearing aids seamlessly adapt to your lifestyle. The Arc AI CROS/BiCROS System enhances audibility and speech understanding, ensuring you can hear and experience the world around you.

## CROS/BiCROS solutions. Now, available with Arc AI.

**Our always-on, always automatic, and always effortless hearing technology is designed so you won't miss out on the conversations and sounds that bring you joy.**



*Audibel's easy-to-use CROS/BiCROS System in rechargeable and non-rechargeable styles delivers our clearest Audibel Sound with 55 million personalized adjustments every hour.*



### CROS

*If you're unable to hear in one ear and have normal hearing in the other ear, you may benefit from a CROS solution. You'll be able to fully experience sound in your weaker-hearing ear by way of your healthy ear.*

#### How it Works

Challenging listening environments now have a solution. Our CROS System comes equipped with a microphone that picks up sounds and voices from the side of your weaker-hearing ear, then transmits them to a hearing aid receiver that's fitted on your healthy ear. This innovative solution ensures you won't miss the important moments and sounds in your life.



### BiCROS

If you're unable to hear in one of your ears and your other ear has hearing loss, you may benefit from a BiCROS solution. This will allow you to experience sound in your ear with little to no hearing and your better ear as well.

#### How it Works

Our BiCROS System sends sound from a microphone, placed in your ear with little to no hearing, to a hearing aid that is worn in your better-hearing ear. This sound can be personalized to your hearing loss, allowing you to hear sounds from both sides.

***Your hearing care professional can help determine what solution will best fit your needs and lifestyle.***